

2024 CATALOG

# **Empowering Healthcare with Advanced Radiation Protection**

#### **Table of Contents**

About INFAB	3
Core Material	4
Revolution Aprons	6
Revolution Apron Sizing Chart	8
Revolution Apron Colors	22
Classic Aprons	26
Classic Apron Sizing Chart	28
Classic Apron Colors	41
Lightning Aprons	52
Smart Track	54
Accessories	56
Scatter Armor Shields	64
C-Arm Drapes	66
Cassette Covers	68
Glasses	70
Gloves	93
Barriers	100
Lockers	110
Apron Racks	112
Curtains	121
Patient Protection	127
Patient Positioning	130
Certification	133
Terms & Conditions	134

### **Engage & Discover**

Join our social channels for the latest in radiation protection and community insights.



# **Linkedin**Network with professionals and see company updates



**Instagram**Follow our visual journey and product highlights.



X (Formerly Twitter) Get the latest news and engage with us.



YouTube
Watch our how-to's and product insights.



**Facebook**Join our community for updates and exclusive content.



Linktree
Scan to find all of our socials in one click.



#### **About Us**

- **Pioneering Advanced Safety:** Designing and manufacturing innovative radiation protection products since 1981.
- American Craftsmanship: Committed to quality with products proudly made in the USA.
- **Technological Breakthroughs:** KIARMOR Bilayer is our lightest, yet most protective lead-free technology to date.
- **Industry Benchmarks:** Consistently elevating the standards for personal radiation protection.

#### **Our Mission & Vision**

- **Healthcare Longevity:** Dedicated to enhancing the welfare of healthcare professionals through superior protection products.
- Innovation Leader: Spearheading the development of lightweight, ergonomic protection that healthcare professionals trust.
- **Global Stewardship:** Fostering a sustainable future with products designed to support healthcare across the globe.

#### **Connect with Our Sales Team**

- Responsive Service: Reach out for personalized assistance Monday - Friday, 6:30 AM - 5 PM (PT).
- **Direct Support Line:** (805) 987-5255 Talk to a specialist who understands your needs.
- Email Inquiry: info@infabcorp.com
   For detailed product information and support.
- Explore Our World: www.INFABCORP.com
  Discover our full range of protective solutions.



# The World's Most Advanced Certified Radiation Protection Core Material

KIARMOR passes all three universally recognized testing procedures:

- ASTM F2547-18
- IEC 61331-1:2014
- DIN 6857-1












# Bilayer: 2 Individual Layers within 1 Homogenous Material

KIARMOR takes advantage of the physics associated with using individual layers of specific materials to provide maximum attenuation by eliminating the effects of fluorescence.





Sb	Antimony	<b>51</b> (Low atomic weight element)
Bi	Bismuth	83 (High atomic weight element)

# **Tested & Certified**

KIARMOR Bi-Layer is tested and certified for attenuation and lead equivalence by the National Physical Laboratory (NPL) against a full range of radiation qualities using the broad beam geometry method.

## **Radiation Protection Core Material**

		Materials	Lead Equivalency	50 kVp	70 kVp	Attenu 90 kVp	aation 110 kVp	120 kVp	150kVp	Core Weight Lb/Ft² (of core in apron)	Core Weight KG/m² (area density)	Front Protection Male Medium Apron Weight 303M-M
	KIARMOR Lead-Free is the lightest-weight protection offered.	Antimony,	.25mm LE	98.7%	96.1%	92.0%	87.8%	85.8%	80.1%	0.522	2.55	3.85 Lbs.
Passes all three univers	Passes all three universally recognized testing procedures:	Bismuth One Homogeneous	.35mm LE	99.6%	98.0%	95.3%	92.7%	91.5%	88.1%	0.729	3.56	5.0 Lbs.
BI-LAYER LEAD-FREE	ASTM F2547-18, IEC 61331- 1:2014, DIN 6857-1	Layer	.5mm LE	99.9%	99.4%	97.8%	96.3%	95.5%	93.0%	1.042	5.09	6.95 Lbs.
4	Our Lightweight Lead		.25mm LE	98.8%	95.1%	90.7%	87.5%	86.2%	82.7%	0.622	3.04	4.35 Lbs
HCHTWEICHT	core material passes all three universally recognized testing procedures: ASTM F2547-18, IEC 61331-1:2014, DIN 6857-1	niversally recognized testing Lightweight Lead ocedures: ASTM F2547-18, IEC	.35mm LE	99.7%	97.6%	94.5%	92.3%	91.3%	88.6%	0.872	4.26	6.15 Lbs
LIGHTWEIGHT LEAD			.5mm LE	99.9%	99.1%	97.4%	96.1%	95.5%	93.0%	1.245	6.08	8.0 Lbs

For over forty years, INFAB has been at the cutting edge of delivering innovative solutions that safeguard against hazardous radiation. Our streamlined selection features two primary core materials: KIARMOR Bilayer Lead-Free and Lightweight Lead, each designed to ensure optimal protection.

Our relentless pursuit of excellence is evident in our research and development efforts, leading to core materials that provide comprehensive attenuation for 50 to 150 kVp. With INFAB's KIARMOR Bilayer, we present the lightest protective material that surpasses all international protection standards, embodying our commitment to safety and quality.

# REVOLUTION

**Key Features** 



**Cool Wear** 



**Abrasion Resistant** 



Flame Resistant



**Water Resistant** 



**Bacteria Resistant** 



**Smart Track** 



Monograms can be added either in "Block" or "Script" to any INFAB Aprons

# How to Measure for Your New Apron

- Take all measurements with a cloth measuring tape.
- Measurements should be taken over your clothing.
- Do not attempt to take your own measurements. You should request the help of an assistant for best accuracy.
- Refer to the Find Your Size diagram for proper location of measurements.
- Some apron patterns are cut for specific genders. In these
  instances, if users of various genders will be sharing
  aprons, plan to consider male size options. If aprons are
  being purchased for separate gender use, plan to select
  the specific gender of the wearer(s) to ensure the best
  protection/fit.
- Once all measurements have been taken, refer to our sizing chart to find your size.
- If it is determined that a user is between sizes, consider the larger size. Garments that are too big may be made smaller; however, garments that are too small cannot be made larger.

# The Tools You Will Need



**Cloth Measuring Tape** 



**Pencil & Paper** 



**Measuring Assistant** 

#### **REVOLUTION | FEMALE SIZING**

#### Revolution Vest Length (\*)

Location	xSmall, Small, Medium, Large, xLarge	2xLarge
Front	22" / 56 cm	24" / 61 cm
Back	26" / 66 cm	28" / 66 cm

#### Revolution Skirt Length (\*)

Location	xSmall, Small, Medium, Large, xLarge, 2xLarge
Wraparound	22" / 56 cm

#### Revolution Front Panel Apron Length (\*)

Location	xSmall, Small, Medium, Large, xLarge, 2xLarge
Front	38" / 97 cm

#### Revolution Full Wrap Apron Length (\*)

Location	xSmall, Small, Medium, Large, xLarge, 2xLarge
Wraparound	38" / 97 cm

#### **Revolution 903 Vest Sizing**


Location	xSmall	Small	Medium	Large	xLarge	2xLarge
Chest (C)	33"-34" / 84-86 cm	35-36" / 89-91 cm	37-38" / 94-97 cm	39-40" / 99-102 cm	41-42" / 104-107 cm	43-44" / 109-112 cm
Length (L)	21" / 53 cm	21.5" / 55 cm	22" / 56 cm	22.5" / 57 cm	23" / 58 cm	23.5" / 60 cm

#### **Revolution 903 Skirt Sizing**

Location	xSmall	Small	Medium	Large	xLarge	2xLarge
Waist (W)	25-27" / 64-67 cm	28-30" / 71-76 cm	31-33" / 79-84 cm	34-36" / 86-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm
Length (L)	22" / 56 cm	22" / 58 cm	22" / 58 cm			

#### **603 Unisex Sizing**

Location	xSmall	Small	Medium	Large	xLarge	2xLarge
Chest (C)	30-33" / 76-84 cm	34-37" / 86-94 cm	38-41" / 97-104 cm	42-45" / 107-114 cm	46-49" / 117-125 cm	50-52" / 127-132 cm
Waist (W)	25-27" / 64-67 cm	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm




#### Female Measurements (Applies to All \*)

Location	xSmall	Small	Medium	Large	xLarge	2xLarge
Chest (C)	30-32" / 76-81 cm	33-36" / 84-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm	43-45" / 109-114 cm	46-48" / 117-122 cm
Waist (W)	25-27" / 64-67 cm	28-30" / 71-76 cm	31-33" / 78-84 cm	34-36" / 86-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm

Length (L): We recommend 1" above the knee.

Contact us for sizes above 2XL



#### Revolution Vest Length (\*)

Location	Small, Medium, Large, xLarge	2xLarge
Front	24" / 61 cm	26" / 66 cm
Back	28" / 71 cm	30" / 71 cm

#### Revolution Skirt Length (\*)

Location	Small, Medium, Large, xLarge, 2xLarge
Wraparound	22" / 56 cm

#### Revolution Front Panel Apron Length (\*)

Location	Small, Medium, Large	xLarge	2xLarge
Front	38" / 97 cm	42" / 107 cm	44" / 112 cm

#### Revolution Full Wrap Apron Length (\*)

Location	Small, Medium, Large, xLarge, 2xLarge
Wraparound	38" / 97 cm

#### **Revolution 903 Vest Sizing**

Location	Small	Medium	Large	xLarge	2xLarge
Chest (C)	37-38" / 94-97 cm	39-40" / 99-102 cm	41-42" / 104-107 cm	43-44" / 109-112 cm	45-46" / 114-117 cm
Length (L)	24.5" / 62 cm	25" / 63.5 cm	25.5" / 65 cm	26" / 66 cm	26.5" / 67 cm

#### **Revolution 903 Skirt Sizing**

Location	Small	Medium	Large	xLarge	2xLarge
Waist (W)	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm
Length (L)	22" / 56 cm	22" / 56 cm	22" / 56 cm	22" / 58 cm	22" / 58 cm

#### Male Measurements (Applies to All \*)

Length (L): We recommend 1" above the knee.

Location	Small	Medium	Large	xLarge	2xLarge
Chest (C)	34-37" / 86-94 cm	38-41" / 97-104 cm	42-45" / 107-114 cm	46-49" / 117-125 cm	50-52" / 127-132 cm
Waist (W)	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm

Contact us for sizes above 2XL

#### 603 Unisex Sizing

Location	Small	Medium	Large	xLarge	2xLarge
Chest (C)	34-37" / 86-94 cm	38-41" / 97-104 cm	42-45" / 107-114 cm	46-49" / 117-125 cm	50-52" / 127-132 cm
Waist (W)	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm

# Revolution Core Material Configurations



Non-Overlapping Sections

## PARTIAL OVERLAP

REVOLUTION Vests and Full Wrap Aprons *DO NOT REQUIRE* a full overlap to achieve the stated protection level shown. Skirt Aprons *REQUIRE* a full overlap to achieve the stated protection level shown.

**MODELS:** 203, B203 MODELS: 103, B103, L103, BL103 **⊢** .25mm LE **⊢** H .35mm LE ⊢ **⊢** .25mm LE **⊢** H .50mm LE ⊢ ط 50mm LE. 50mm LE 70mm LE 1.0mm LE 50mm LE 70mm LE 1.0mm LE .50mm LE .50mm LE All non-overlapping skirt Overlapping panels offer extra protection. Standard Standard panels are .25mm LE Configuration Configuration

# Revolution Core Material Configurations





## **NO OVERLAP**


REVOLUTION Front Wrap Aprons achieve the stated protection level shown.

## **FULL OVERLAP**

REVOLUTION 703 and 803 Aprons *REQUIRE* a full overlap to achieve the stated protection level shown.


**MODELS:** 303, B303, 403, B403, 503, B503, 603, B603

MODELS: 703, 803



S0mm LE

.50mm LE



The 603 and B603 back panels are .25mm LE

.50mm LE Standard Configuration

.50mm LE Standard Configuration

All non-overlapping panels are .25m LE

#### **REVOLUTION**













#### **Revolution Lumbar Vest & Skirt**

- Built-in "Black Belt" in vest and skirt for better lumbar support
- Partial overlap to increase protection
- Features interior pocket on both vest and skirt
- "Back Flap" for better weight balance
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack Track enabled

Model	Description
L103	Revolution Lumbar Vest and Skirt
L103V	Revolution Lumbar Vest
L103S	Revolution Lumbar Skirt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-2XL Female XS-2XL Length alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	BL103 Classic Nylon Colors without Cool Wear
Monogram	Block or <i>Script</i>















Model	Description
103	Revolution Vest and Skirt
103V	Revolution Vest
103S	Revolution Skirt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-2XL Female XS-2XL Length alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B103 Classic Nylon Colors without Cool Wear
Monogram	Block or Script

Color Shown:

# 103 **Revolution Vest & Skirt**

- Features interior pocket on both vest and skirt
- "Back Flap" for better weight balance
- Removable Black Belt® in skirt for better lumbar support
- Partial overlap to increase protection
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

#### **REVOLUTION**












# 203

## **Revolution Full Wrap Black Belt**

- Partial overlap to increase protection
- Hook and loop closure side slits for increased mobility
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- Snap on, fully adjustable Black Belt® for lumbar support
- SmartTrack enabled

Model	Description
203	Revolution Full Wrap Black Belt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35mm LE, .50mm LE
Sizes	Male S-2XL Female XS-2XL Length alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B203 Classic Nylon Colors without Cool Wear
Monogram	Block or Script











Model	Description
303	Revolution Front
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-2XL Female XS-2XL Length alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B303 Classic Nylon Colors without Cool Wear
Monogram	Block or <i>Soript</i>

# **303**

#### **Revolution Front**

- Criss-cross back straps to provide weight relief and lumbar supportCool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

#### **REVOLUTION**










Color Shown: 513 Tide Red - Black Block

# 403

#### **Revolution Front Elastic Back**

- Criss-cross elastic back straps to provide weight relief and lumbar support
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled


Model	Description
403	Revolution Front Elastic Back
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-2XL Female XS-2XL Length alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B403 Classic Nylon Colors without Cool Wear
Monogram	Block or Script







Color Shown: 504 Pink - Black Block






# **503**


#### **Revolution Front Black Belt**

- Cool Wear moisture-wicking interior fabric
- 503 comes with a snap on, fully adjustable Black Belt®
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

#### **REVOLUTION**









Color Shown: 602 Black - Fuchsia Block

# 603

## **Revolution Full Wrap Elastic Back**

- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

Model	Description
603	Revolution Full Wrap Elastic Back
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Unisex XS-2XL Length Alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B603 Classic Nylon Colors without Cool Wear
Monogram	Block or <i>Soript</i>







Color Shown: 510 Purple - Black Block






# **703**

#### Revolution Full Overlap Lumbar Vest & Skirt

- Full overlap
- Built-in Black Belt® in vest and skirt for better lumbar support
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

#### **REVOLUTION**









**Color Shown:** 515 Orange - Black Block

# 803

## Revolution Black Belt® Wrap

- Full overlap
- Hook and loop closure side slits for increased mobility
- Cool Wear moisture-wicking interior fabric
  Professionally designed for an exceptional fit and comfort
- Snap on, fully adjustable Black Belt®
- SmartTrack enabled

Model	Description
803	Revolution Black Belt Wrap
Colors	500 Series, 600 Series, 700 Series, 800 Series, 900 Series, Remix
Core	Lightweight Lead, KIARMOR Lead-Free
Protection	.50mm LE
Sizes	Female XS-2XL Male S-2XL Length Alterations ±2"or ±4"
Fabric	Revolution Fabric with Cool Wear
Monogram	Block or <i>Soript</i>







Color Shown: 512 Ducks Green - Black Block





## 903

#### **Revolution Reverse Vest & Skirt**

- Criss-cross elastic back straps to provide weight relief and lumbar support
- Adjustable shoulder strap ensures a precise fit
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

# **Revolution Colors**

## Color

**700 SERIES** 



## Color + Block

**500 SERIES** 



Color Shown: 506



511 White -Black Block

## **Black + Color**

#### **600 SERIES**



Color Shown: 606



611 Black -

White Block

## **Color + White**

#### **800 SERIES**



Color Shown: 806





# **Revolution Colors**

## White + Color

900 SERIES



Color Shown: 906



White -Black Block

# Color + Color

**REMIX SERIES** 





# **Revolution Base**



BL103 Vest & Skirt Lumbar Belt



**B103** Vest & Skirt



**B203**Full Wrap Apron Lumbar Belt



**B303**Front Panel Apron



**B403**Front Panel Apron - Elastic Back



**B503**Front Panel Apron Lumbar Belt



**B603**Full Wrap Apron Elastic Back

Revolution Base offers all the comfort and protection of our premium Revolution line, but with our classic nylon colors and patterns. Revolution Base does not include Revolution Fabric, Color-Blocking or Cool-Wear.

Size Female XS-2XL

Male S-2XL

Fabric Nylon

Protection .25mm LE, .35mm LE, and .50mm LE (Options Vary per Model)

Core Lightweight Lead or KIARMOR Lead-Free
Color Classic color & patterns options only

# **Classic Aprons**

**Key Features** 



**Water Resistant** 



**Abrasion Resistant** 



**Flame Resistant** 



**Smart Track** 



**Bacteria Resistant** 



**Custom Patterns** 



Monograms can be added either in "Block" or "8cript" to any INFAB Aprons

# How to Measure for Your New Apron

- Take all measurements with a cloth measuring tape.
- Measurements should be taken over your clothing.
- Do not attempt to take your own measurements. You should request the help of an assistant for best accuracy.
- Refer to the Find Your Size diagram for proper location of measurements.
- Some apron patterns are cut for specific genders. In these
  instances, if users of various genders will be sharing
  aprons, plan to consider male size options. If aprons are
  being purchased for separate gender use, plan to select
  the specific gender of the wearer(s) to ensure the best
  protection/fit.
- Once all measurements have been taken, refer to our sizing chart to find your size.
- If it is determined that a user is between sizes, consider the larger size. Garments that are too big may be made smaller; however, garments that are too small cannot be made larger.

# The Tools You Will Need



**Cloth Measuring Tape** 



**Pencil & Paper** 



**Measuring Assistant** 

# Find Your Size

#### **HELPFUL HINTS**

Chest (C): Measure widest part 1" below armhole.

Waist (W): Measure narrowest part 1" above navel.

Hips (H): Measure widest part of hips.

Length (L): We recommend 1" above the knee.

#### Classic Front Panel Apron Length (\*)

- 1	Location xSmall, Small, Medium, Large, xLarge, 2			
Front		38" / 97 cm		

#### Classic Full Wrap Apron Length (\*)

Location	xSmall, Small, Medium, Large, xLarge, 2xLarge
Wraparound	38" / 97 cm

Measure from high point shoulder to about 1" above the knee. This will be your one piece apron length. If you see this icon, it means that the apron has its own unique measurements.

**SPECIAL** 

Measure from high point shoulder to desired skirt belt line. Then add 4" to find your vest length.

Measure from your desired skirt belt line to about 1" above the knee. This will be your skirt length.

#### Reverse & Elastic Vest Sizing (\*)

Location	xSmall	Small	Medium	Large	xLarge	2xLarge
Length (L)	24" / 61 cm	24" / 61 cm	26" / 66 cm	27" / 69 cm	28"/71 cm	29" / 74 cm

#### Classic Vest Sizing (\*)

Location	xSmall	Small	Medium Large		xLarge	2xLarge
Length (L)	20" / 51 cm	22" / 56 cm	24" / 61 cm	25" / 64 cm	26" / 66 cm	27" / 69 cm

#### Classic Skirt Sizing (\*)

Location	xSmall	Small	Medium Large		xLarge	2xLarge
Length (L)	18" / 47 cm	20" / 51 cm	22" / 56 cm	23" / 58 cm	24" / 61 cm	25" / 64 cm

Available Length
Alterations in increments of ±1"


#### Female Measurements (Applies to All \*)

(C)

(W)

Location	on xSmall Small Medium Large		xLarge	2xLarge		
Chest (C)	30-32" / 76-81 cm	33-36" / 84-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm	43-45" / 109-114 cm	46-48" / 117-122 cm
Waist (W)	25-27" / 64-67 cm	28-30" / 71-76 cm	31-33" / 78-84 cm	34-36" / 86-91 cm	37-39" / 94-99 cm	40-42" / 102-107 cm

Contact us for sizes above 2XL



# **Find Your Size**

#### **HELPFUL HINTS**

Chest (C): Measure widest part 1" below armhole.

Waist (W): Measure narrowest part 1" above navel.

Hips (H): Measure widest part of hips.

Length (L): We recommend 1" above the knee.

#### Classic Front Panel Apron Length (\*)

Location	Small, Medium, Large	xLarge	2xLarge
Front	38" / 97 cm	42"/107 cm	44"/112 cm

#### Classic Full Wrap Apron Length (\*)

Location	Small, Medium, Large, xLarge, 2xLarge
Wraparound	38" / 97 cm

#### Reverse & Elastic Vest Sizing (\*)

Location	Small	Medium	Large	xLarge	2xLarge
Length (L)	25" / 64 cm	27" / 69 cm	28" / 71 cm	30" / 76 cm	30" / 76 cm

#### Classic Vest Sizing (\*)

apron length.

Location	Small	Medium	Large	xLarge	2xLarge
Length (L)	22" / 56 cm	24" / 61 cm	25" / 64 cm	26" / 66 cm	28" / 71 cm

#### Classic Skirt Sizing (\*)

Location	Small	Medium	Large	xLarge	2xLarge
Length (L)	20" / 51 cm	22" / 56 cm	23" / 58 cm	24" / 61 cm	25" / 64 cm

#### Male Measurements (Applies to All \*)

Location	Small	Medium	Large	xLarge	2xLarge
Chest (C)	34-37" / 86-94 cm	38-41" / 97-104 cm	42-45" / 107-114 cm	46-49" / 117-125 cm	50-52" / 127-132 cm
Waist (W)	28-31" / 71-79 cm	32-35" / 81-89 cm	36-39" / 91-99 cm	40-43" / 102-109 cm	44-46" / 112-117 cm

Contact us for sizes above 2XL

# Classic Core Material Configurations




Non-Overlapping Sections

# **FULL OVERLAP**


CLASSIC Full Wrap Vest and Skirt Aprons REQUIRE a full overlap to achieve the stated protection level shown.

#### **CLASSIC:** VAS

#### **CLASSIC:** LSW, FCW



.50mm LE Standard Configuration



.50mm LE Standard Configuration

# Classic Core Material Configurations






## **NO OVERLAP**

CLASSIC Front Wrap Apron's and Vest's achieve the stated protection level shown.

Only the CLASSIC Skirt *REQUIRE* a full overlap to achieve the stated protection level shown.

CLASSIC: AWC, EBA, LSF, MWA, WBT, SPL, SDO, SDA, TFAM, DTL-TECH, DTL-PANO, DTL-ADULT, DTL-ADULT-TC, DTL-CHILD, DTL-CHILD-TC • EVAS, RVAS



#### **VAS** Vest and Skirt









The Classic Vest and Skirt apron provides full overlap protection and maximum weight distribution between the shoulders and hips, eliminating stress on the upper and lower back. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29)

FABRIC: Nylon (with a double welt pocket)

**PROTECTION: .50mm LE** 

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

#### **RVAS** Reverse Vest & Skirt











The Reverse Vest and Skirt apron provides maximum weight distribution. Criss-cross elastic back straps provide weight relief and lumbar support, while the adjustable shoulder strap ensures a precise fit. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### **EVAS** Elastic Vest & Skirt









The Elastic Vest and Skirt apron provides maximum weight distribution. Criss-cross elastic back straps provide weight relief and lumbar support, while the adjustable back strap ensures a precise fit. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### **LSW** Full Overlap with Lumbar Belt









The Black Belt Full Wrap apron offers maximum protection with full overlap, featuring a snap-on, fully adjustable lumbar Black Belt®. This apron combines comfort with superior lumbar support, which reduces fatigue and upper back stress during lengthy procedures. Built-in slits on each side with hook and loop closures provide increased mobility. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29)

FABRIC: Nylon (with a double welt pocket)

**PROTECTION:** .50mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

# MWA Special Wrap







The Special Wrap apron is uniquely designed to offer full protection with an adjustable shoulder strap and front belt buckle closure. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### **FCW** Full Coat Wrap









The Full Coat Wrap offers maximum front and back full overlap protection. Built-in slits on each side with hook and loop closures provide increased mobility. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29)

FABRIC: Nylon (with a double welt pocket)

**PROTECTION:** .50mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

#### **AWC** Weight Reliever







The Weight Reliever front protection apron features side closure straps at the waist which relieve weight from the shoulders and adjust for the most comfortable fit. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29) **FABRIC:** Nylon (with a double welt pocket) PROTECTION: .25mm LE, .35mm LE, and .50mm LE CORE: Lightweight Lead, KIARMOR Lead-Free PRINT: Solid and pattern options (Pages 41-51)

#### **EBA Elastic Back**







The Elastic Back front protection apron features elastic side closure straps providing weight relief, comfort, and flexibility. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29) FABRIC: Nylon (with a double welt pocket) **PROTECTION:** .25mm LE, .35mm LE, and .50mm LE CORE: Lightweight Lead, KIARMOR Lead-Free PRINT: Solid and pattern options (Pages 41-51)

#### **Black Belt Front**





The Black Belt Front Protection apron provides frontside body protection while offering snap-on, fully adjustable Black Belt® lumbar support. The built-in back strap provides better weight distribution between the shoulders and waist. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29) **FABRIC:** Nylon (with a double welt pocket) PROTECTION: .25mm LE, .35mm LE, and .50mm LE CORE: Lightweight Lead, KIARMOR Lead-Free **PRINT:** Solid and pattern options (Pages 41-51)

#### **WBT Wide Belt**









The Wide Belt apron offers frontside protection with a belt and buckle closure, adjustable for a comfortable fit. SmartTrack enabled.

SIZE: xSmall - 2xLarge (Pages 28-29) **FABRIC:** Nylon (with a double welt pocket) PROTECTION: .25mm LE, .35mm LE, and .50mm LE CORE: Lightweight Lead, KIARMOR Lead-Free **PRINT:** Solid and pattern options (Pages 41-51)

#### **SPL** Budget Saver







The Budget Saver front protection apron features criss-cross tie straps for a secure, snug fit and comfortable shoulder pads. SmartTrack enabled.

**SIZE:** xSmall - 2xLarge (Pages 28-29)

FABRIC: Nylon

PROTECTION: .25mm LE, .35mm LE, and .50mm LE
CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

#### **SDO** Surgery Drop Off







The Surgery Drop Off front protection apron is designed to be removed without breaking the sterile field. Additional handle straps located on top of shoulders for easy removal or drop off. Attaches in the front with hook and loop closures. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### **SDA** Surgery Drop Away







The Surgery Drop Away front protection apron features hook and loop closures at the top of each shoulder for easy removal without breaking the sterile field. The back strap provides a secure fit. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### TFA Tri-Flap







The Tri-Flap apron features a split front panel to provide protection and comfort while seated. Criss-cross straps attach in front with a hook and loop closure. SmartTrack enabled.

SIZE:xSmall - 2xLarge (Pages 28-29)FABRIC:Nylon (with a double welt pocket)PROTECTION:.25mm LE, .35mm LE, and .50mm LECORE:Lightweight Lead, KIARMOR Lead-FreePRINT:Solid and pattern options (Pages 41-51)

#### **DTL-TECH** Dental Technician Apron

#### **DTL-PANO** Pano Cape Dental Apron





The Dental Technician Apron front protection apron features criss-cross tie straps for a secure, snug fit and comfortable shoulder pads. SmartTrack enabled.

**SIZE:** 24" W x 36" L

**FABRIC:** Nylon **PROTECTION:** .35mm LE

**CORE:** Lightweight Lead, KIARMOR Lead-Free **PRINT:** Solid and pattern options (Pages 41-51)



The Pano Cap Dental Apron drapes over your patient's shoulders. The inner liner of this cape uses a non-slip material and fastens in the back with a hook and loop closure. SmartTrack enabled.

**SIZE:** 23" W x 25" L

**FABRIC:** Nylon **PROTECTION:** .35mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

#### **DTL-ADULT** Adult Dental Apron



The Adult Dental Apron provides front protection for the patient, using a non-slip material to keep it firmly in position. An attached thyroid collar option is available in this style. SmartTrack enabled.

**SIZE:** 24" W x27" L

**FABRIC:** Nylon **PROTECTION:** .35mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

#### DTL-CHILD Child Dental Apron



The Child Dental Apron provides front protection for the patient, using a non-slip material to keep it firmly in position. An attached thyroid collar option is available in this style. SmartTrack enabled.

**SIZE:** 18" W x 24" L

**FABRIC:** Nylon **PROTECTION:** .35mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free
PRINT: Solid and pattern options (Pages 41-51)

## **Apron Colors**

We use high quality nylon and polyester fabrics for all of our products. A polyurethane coating seals the cloth from fluid penetration and makes it easier to clean.

Color-fast dyes won't transfer to your clothing. Colors and patterns may be added and discontinued at any time.



#### **Nature**





SH-20 New Flames



SH-32 Orange Hibiscus



180 Orange Hawaii



LT222-37 Perfect Petals



LT221-16 Lovely Lemons



SH-49 Tropical Avocado



LT221-22 Spring Garden



MIR272-53 Mantis Blue



217: Mardi Gras



SH-19 Pink Paisley




**SH-28** 80's Palm



SH-13 Hawaiian Flowers



SH-24 Cactus



208 Blue Wash



207 Blue Flame



SH-54 Lightning




SCP244-85 Teal Splash



LT222-35 Fancy Foliage

SH-31 Sugar Skulls



All Patterns Available on Classic & Revolution Base Aprons

140 Splash



SH-22 Fitness



SH-17 Race Car

### **Sports**

185 White Dots on Pink



265 Chicklets Purple

SH-21 Volleyball



209 Chicklets

SH-9: Sports


# Camo 201 Digital Camo 268 Digital Camo Black 120 Camouflage SH-18 Shape Camo MIR215-51 Blue Camo MIR215-55 Navy Camo 184 Pink Camouflage SH-6 Color Camo TRTB-NC65 Timber Camo SH-27 Orange Blocks **SH-42** 80's **Abstract** SH-7 Metallic Splatter SH-2 Blocks

LT222-36 Marble Geo

LT221-13 Gold Geo

176 Rainbow Wave

SH-23 Art Deco



### Wildlife



SH-52 Mixed Dogs



SH-40 Frenchies



SH-44 80's Abstract Frenchies



SH-37 Dachshund



SCP304-09 Dog Paw



167 Cheetah



CAS009-33 Wildcat Pink



214 Animal Print Pink



212 Animal Print Blue



SH-10 Cats



266 Tiger Skin



263 Chocolate Giraffe



SH-50 Colorful Cats



SH-26 Peacock



SH-55 Ocean Life



SH-30 Sharks



SH-8 Fish Bones



262 Black Butterflies



SH-15 Flamingo

SH-12 Unicorn



216 Dragons

We use high quality nylon and polyester fabrics for all of our products. A polyurethane coating seals the cloth from fluid penetration and makes it easier to clean.

SH-36 Colored Fossils

SH-14 Dino Fossil

Color-fast dyes won't transfer to your clothing. Colors and patterns may be added and discontinued at any time.

SH-48 Black Unicorn

For the most up to date list of apron colors, visit www.infabcorp.com/colors

#### Leaving Soon

**Check Availability** infabcorp.com/colors



172 Puppies



SCP145-95 Red Granite



TFX767-37 Baby Leo





219 Pink Paws



218 Mamba



SCP088-43 Purple Flame





264: Sheep in Clouds Pink



SCP305-95 Dalmatians



TXF779-44 Mystery Splash





TFX755-31 Pink Plaid



MIR260-03 Candy Circles



MIR239-53 Happy Talk



SH-5 Teeth




MIR272-33 Mantis Pink



TFX755-45 Purple Plaid



SCP258-54 Vintage Postcard



For the most up to date list of apron colors, visit www.infabcorp.com/colors

## **CLASSIC Ripstop**

Ripstop fabrics are woven fabrics, often made of nylon, using a special reinforcing technique that makes them resistant to tearing and ripping. A polyurethane coating seals the cloth from fluid penetration and makes it easier to clean. Color-fast dyes won't transfer to your clothing.





Model: B103 Print Shown: 194 Burgundy Diamond Ripstop

### **Trim Colors**

INFAB offers a variety of colors for your apron. Color-fast dyes won't transfer to your clothing. Colors and patterns may be added and discontinued at any time. For the most up to date list of apron colors, visit www.infabcorp.com/colors







## Monogram

INFAB offers customized monograms in Block or Script styles to personalize your radiation protection aprons.

#### **Block**

AaBbCcDdEeFfGg HhliJjKkLIMmNnOo PpQqRrSsTtUuVv WwXxYyZz 1234567890

### Script

AaBb Cc Dd Ee Ff Gg Hh Ii Jj Kk Ll Mm Nn Oo Pp Qq Rr Ss Tt Uu Vv Ww Xx Yy Zz 1234567890

# LIGHTMING

### **Quick Ship Aprons**







**Smart Track** 



**Abrasion Resistant** 



**Flame Resistant** 



**Water Resistant** 



**Bacteria Resistant** 

### Available in Six Colors



114: Black



108: Royal Blue



205: Charcoal



119: Purple



105:

Red

Teal

117:

#### **Lightning Vest & Skirt**



Model: UVS Color Shown: Purple

The Lightning Vest and Skirt apron provides partial overlap protection and maximum weight distribution between the shoulders and hips, eliminating stress on the upper and lower back. SmartTrack enabled.

SIZE: Small - 2xLarge

**FABRIC:** Nylon (with a double welt pocket)

PROTECTION: .35mm LE, and .50mm LE

CORE: KIARMOR Lead-Free, Lightweight Lead, **COLORS:** Black, Blue, Charcoal, Purple, Red, Teal

Size	Chest	Waist	Vest Front	Skirt
S	34-37" / 86-94 cm	28-31" / 71-79 cm	24" / 61 cm	20" / 51 cm
М	38-41" / 97-104 cm	32-35" / 81-89 cm	24" / 61 cm	22" / 56 cm
L	42-45" / 107-114 cm	36-39" / 91-99 cm	24" / 61 cm	23" / 58 cm
XL	46-49" / 117-125 cm	40-43" / 102-109 cm	24" / 61 cm	24" / 61 cm
2XL	50-52" / 127-132 cm	44-46" / 112-117 cm	24" / 61 cm	25" / 64 cm

#### **UFW** Lightning Full Wrap

The Lightning Full Wrap offers maximum front and back partial overlap protection. Built-in slits on each side with hook and loop closures provide increased mobility. SmartTrack enabled.

**SIZE:** Small - 2xLarge

**FABRIC:** Nylon (with a double welt pocket)

**PROTECTION:** .35mm LE, and .50mm LE

**CORE:** KIARMOR Lead-Free, Lightweight Lead, **COLORS:** Black, Blue, Charcoal, Purple, Red, Teal

Size	Chest	Waist	Length
S	34-37" / 86-94 cm	28-31" / 71-79 cm	38" / 97 cm
М	38-41" / 97-104 cm	32-35" / 81-89 cm	38" / 97 cm
L	42-45" / 107-114 cm	36-39" / 91-99 cm	38" / 97 cm
XL	46-49" / 117-125 cm	40-43" / 102-109 cm	38" / 97 cm
2XL	50-52" / 127-132 cm	44-46" / 112-117 cm	38" / 97 cm

#### **UFP** Lightning Front Protection





The Lightning Front Protection apron features side closure straps at the waist which relieve weight from the shoulders and adjust for the most comfortable fit. SmartTrack enabled.

SIZE: Small - 2xLarge

**FABRIC:** Nylon (with a double welt pocket)

**PROTECTION:** .35mm LE, and .50mm LE

CORE: KIARMOR Lead-Free, Lightweight Lead,
COLORS: Black, Blue, Charcoal, Purple, Red, Teal

Size	Chest	Waist	Length
S	34-37" / 86-94 cm	28-31" / 71-79 cm	38" / 97 cm
М	38-41" / 97-104 cm	32-35" / 81-89 cm	38" / 97 cm
L	42-45" / 107-114 cm	36-39" / 91-99 cm	38" / 97 cm
XL	46-49" / 117-125 cm	40-43" / 102-109 cm	38" / 97 cm
2XL	50-52" / 127-132 cm	44-46" / 112-117 cm	38" / 97 cm



# Track & Document Lead Apron Inspections •——

Manage all of your lead aprons in one place with Smart Track, the completely FREE and secure web-based garment tracking system.



#### **HOW IT WORKS**

- 1. Download the free app.
- 2. Create your Smart Track account.
- 3. Add your garments.
- 4. Document your inspections.
- 5. Set inspection reminders.
- 6. Export and download reports.

#### **KEY FEATURES**

- No Software & No Registration
- · Garment Dashboard
- · Custom Reports
- Bulk Uploads & Editing
- Filtered Search
- + More!







#### TRACK ANY GARMENT

Use your existing garment IDs or add Smart Track Retrofit Tags to any apron to get started. Smart Track Retrofit Tags feature unique IDs that can be scanned with the Smart Track App or any 2D scanner.









SMART TRACK Scanner





SMART TRACK 2D bar code tags with permanent securitystyle backing

#### **TESTIMONIALS**

"Smart Track makes it simple to manage my suite of lead aprons. We track over 100 aprons - it saves time and headaches."

Rich Norman – Lead Tech Los Robles Hospital, Thousand Oaks, CA

"Tracking lead aprons and making sure that they have been inspected is a challenge for any large radiology facility. Of the many apron-tracking systems I reviewed, I knew that the "SMART TRACK" apron inventory system would work for our facility and meet new regulatory standards set by "The Joint Commission"...The system is user friendly. All I know is that anyone who has the burden of tracking aprons needs to get "Smart" about tracking their apron inventory".

Roseann Zeldin, RT, (R) (M), CNMT – Assistant Radiation Safety Officer Department of Radiology & Imaging Hospital for Special Surgery New York City, New York

"The SMART TRACK Apron Inventory system has revolutionized the way lead apron quality control testing and inventory management is performed at our facility. It is a superior, comprehensive, and user-friendly product. Individual bar make tracking aprons quick and easy. It is the best system I have had the opportunity to work with."

Doug Beardsall MRT (R) CTIC – Senior Technologist QC Royal Victoria Regional Health Center Barrie, Ontario



#### YouTube Videos

Access our helpful videos on how to use the Smart Track App



#### **REV-TC**

#### **Revolution Thyroid Collar**

**SIZES:** One size fits most

**COLORS:** 500 Series, 600 Series, 700 Series, 800 Series, 900 Series, Remix

**CORE:** Lightweight Lead, KIARMOR Lead-Free

**CLOSURE:** Hook and Loop Closure, Buckle

· Large coverage area

• Exclusive Revolution outer fabric

- · Cool Wear moisture-wicking liner
- Optional monogram
- Optional KTC attachment



#### **SVT**

#### **Visor Style Thyroid Collar**

SIZES: X-Small, Standard, Medium

COLORS: A

**CORE:** Lightweight Lead, KIARMOR Lead-Free

**CLOSURE:** Hook and Loop Closure, Buckle

- · Comfortable fit
- · Large coverage area
- Optional monogram
- · Optional KTC attachment



#### TCV

#### **Thyroid Collar with Binding**

SIZES: X-Small, Standard, Medium

**COLORS:** Revolution 700 series & All Classic Nylon Colors

CORE: Lightweight Lead, KIARMOR Lead-Free
CLOSURE: Hook and Loop Closure, Buckle or Magnet

· Comfortable fit

- Large protection area
- · Optional monogram
- · Optional KTC attachment



#### MTC

#### **Magnetic Thyroid Collar**

SIZES: X-Small, Standard, Medium

COLORS: All

**CORE:** Lightweight Lead, KIARMOR Lead-Free

CLOSURE: Magnet

- Comfortable fit
- · Large coverage area
- Optional monogram
- Optional KTC attachment







### NST Thyroid Collar Without Binding

SIZES: X-Small, Standard, Medium

**COLORS:** Revolution 700 Series & All Classic Nylon Colors

CORE: Lightweight Lead, KIARMOR Lead-Free
CLOSURE: Hook and Loop Closure, Buckle, Magnet

· Comfortable fit

- Large protection area
- · Optional monogram
- Optional KTC attachment





KTC Attachment Options



**KTC** Buckle



**KTC** Key Ring



**KTC** No Buckle



KTC attachment is now securely positioned on the shoulder of all Revolution aprons.


300102 Thyroid Collar Cover



- Keeps thyroid collars sanitary and clean
- · Attaches easily to NST-S1 & TCV-S1 collars
- · Easy to wash or replace when soiled
- Revolution 700 Series & All Classic nylon colors
- Can be monogrammed for personalization and easy identification



#### NSM Non-Slip Material for Skirts



- Prevents your INFAB skirt from slipping down the waist and causing discomfort
- Can be added to any INFAB skirt by adding part number NSM



#### **Shin Guards**

- Reduce radiation exposure to the shins
- Hook and loop closures for easy adjustments

**COLORS:** Classic nylon colors **SIZES:** Small, Medium, Large

**CORE:** Lightweight Lead, KIARMOR Lead-Free

#### **REV-BAG** Revolution Apron Bag



- Revolution fabric
- Plastic bottom to prevent sagging
- Opening on top of bag
- · Zipper across the top of bag



#### 909821 Apron Bag



#### **REV-SLV-1** Revolution Attached Thyroid Collar Sleeves



- Revolution outer fabric
- Additional protection for chest
- Wear under any apron

**COLORS:** Revolution 700 Series

SIZES: Small, Medium, Large, X-Large, XX-Large

CORE: Lightweight Lead, KIARMOR Lead-Free



#### **REV-SLV-2** Revolution Sleeve Individual



- Revolution outer fabric
- Hook and loop closure elastic strap for easy on/off convenience
- One per order

**COLORS:** Revolution 700 Series

CORE: Lightweight Lead, KIARMOR Lead-Free





### **SLV** Optional Left or Right Sleeve

- · Added protection of the arm and shoulder
- .50mm LE
- · Available on most INFAB apron models
- Optional removable snap on, permanent sewn on or hook and loop closure attached configuration

COLORS: Revolution 700 Series & All Classic Nylon Colors

ATTACHED: Sewn on, snaps or hook and loop closure CORE: Lightweight Lead, KIARMOR Lead-Free



#### 257890

#### **Suspenders**

- Fully adjustable suspenders with a handy snap on style.
- · One size fits all.



#### The Original Black Belt®

The Black Belt® is a 6" wide orthopedic elastic belt that supports your lower back while removing the weight of the apron from your shoulders. Optional KTC Buckle Attachment and/or side snaps.

**COLORS:** Black

SIZES: X-Small, Small, Medium, Large, X-Large, XX-Large

257879 6" Black Belt® - X-Small 257880 6" Black Belt® - Small 257881 6" Black Belt® - Medium 257882 6" Black Belt® - Large 257883 6" Black Belt® - X-Large 257884 6" Black Belt® - XX-Large



#### **APW** | Axillary Protection Wings

Increase the protection to your axillary region without sacrificing any mobility or adding weight to your apron.

**COLORS:** Revolution 700 Series & All Classic Nylon Colors

CORE: Lightweight Lead, KIARMOR Lead-Free



#### **Revolution Elastic Back Thinking Cap**

• Revolution outer fabric • Cool Wear moisture-wicking fabric

• Scatter radiation protection • Elastic-back for easy on/off convenience

**SIZE:** Small, Medium, Large, X-Large

COLORS: 700 Series Revolution Colors

PROTECTION: .25mm LE, .35mm LE, .50mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free





#### **Classic Thinking Cap with Tie Straps**

Classic nylon colors
 Scatter radiation protection

Head and skull radiation protection
 Cool Wear moisture-wicking fabric

SIZE: Small, Medium, Large, X-Large, XX-Large

COLORS: Classic Nylon Colors

PROTECTION: .25mm LE, .35mm LE, .50mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free





## 400025 Revolution Hook and Loop Closure Back Thinking Cap

Revolution outer fabric

• Cool Wear moisture-wicking fabric

Scatter radiation protection
 Hook & loop back closure for customizable fit

SIZE: One Size Fits All

COLORS: 700 Series Revolution Colors

PROTECTION: .25mm LE, .35mm LE, .50mm LE

CORE: Lightweight Lead, KIARMOR Lead-Free





#### DC-25 Lead-Free Disposable Thinking Cap

· Quantity 20/Case

• .25mm LE

· Adjusts to fit all head sizes

For use by all personnel during any interventional procedure

• 100% lead free for easy disposal

SIZE: One Size Fits All
COLORS: Scatter Armor Blue

**PROTECTION:** .25mm LE

**CORE:** KIARMOR Lead-Free







#### PARTNERING TO KEEP YOUR APRONS CLEAN & SAFE

INFAB is proud to be a global cleaning solutions partner with Clorox.



#### **Excellent product aesthetics**

No harsh chemical odor or fumes

No fragrance added

Nonflammable, unlike quaternary ammonium-alcohol-based wipes

No volatile organic compounds (VOCs)

#### **Broad healthcare surface compatibility**

Recommended for plastics and laminated surfaces, whereas isopropyl alcohol will attack some forms of plastic, rubber and coatings\*

Surface materials: All INFAB aprons, stainless steel, vinyl, glass, baked enamel, chrome, glazed tile, glazed ceramic, glazed porcelain, porcelain enamel, laminated surfaces, plastic If light residue appears, it can easily be removed with a clean, damp cloth



<sup>\*</sup> Based on 30-second non-TB, non-Staphylococcus epidermidis bacteria and virus contact times for general disinfecting on federal master labels as of 2/14/13.



#### 683447

#### Revolution Scrubbles Apron Cleaner

Scrubbles apron cleaner safely cleans your aprons of barium, blood and other common radiology related stains. Now with a clean citrus scent and a stronger more efficient formula.

SIZE 32 oz

Also Available 683448

Case of 12 32 oz bottles with Spray Caps



#### CLX28

#### **Clorox Spray**

Recommended for: Multipurpose use on hard and soft surfaces

SIZE 32 oz



#### CLX24

#### Clorox Multipurpose Wipes

Recommended for: Everyday cleaning and disinfecting, daily patient room cleaning, hard, high-touch surfaces

**SIZE** 6.75" x 9" **OTY** 95 ct



#### CLX26

#### Clorox Terminal Wipes

Recommended for: Terminal cleaning and disinfecting of hard surfaces

**SIZE** 12" x 11" **QTY** 185 ct.



#### CLX27

## Clorox Terminal Wipes Refill

Recommended for: Terminal cleaning and disinfecting of hard surfaces

**SIZE** 12" x 11" **QTY** 185 ct.



#### **Scatter protection**

Shields provide scatter radiation protection for radial and many other types of procedures. KIARMOR Bi-Layer is the most protective core material against scatter radiation.

#### **How it works**

Sterile Scatter Armor Shields are placed on top of the patient drape to block scatter radiation at the source, creating a scatter-free zone, reducing harmful radiation exposure to the doctor or clinician.

#### One-time use

Scatter Armor Shields are lead-free and disposable.

#### Sterile

Packaged individually in sterile packets.

#### **Flexible**

Scatter Armor Shields easily mold to the patients body for ease of use and maximum protection

### **Scatter Radiation Testing According to ASTM Test Method F2547-18.**

Sample Information	50 kV	70 kV	90 kV	110 kV	120 kV	150 kV	
Scatter Armor Shield	98.7%	96.1%W	92.0%	87.8%	85.8%	80.1%	

#### **Direct Beam Attenuation**

Sample Information	60 kV	80 kV	100 kV	120 kV
Scatter Armor Shield	0.967	0.911	0.846	0.788

#### **Absorbent**


Scatter Armor Shields utilize the highest absorbency covering material in the industry, critical in heavy fluid situations.

Absorbent Rate	1.3 sec	
Weight	2.36 osy / 80 gsm	
Absorbent Capacity	686% 16.2osy / 549 gsm	

GSM - Grams per square meter, OSY - Ounce per square yard


#### Scatter Armor Multipurpose Shield

Transradial Access: Left and Right Radial, Shunt, Biliary Drainage, Spinal, Fibroid Embolization




#### Scatter Armor Multipurpose Shield with Fenestration

Needle Biopsy or Puncture Procedures.




#### Scatter Armor Angiography Shield

Transfemoral Access: Angiography, Coronary Catheterization



#### **Scatter Armor EP Shield**

Subclavian Access: Bi-Ventricular Pacemaker Implants, ICS procedures, and Bi-Ventricular Pacing.




#### Scatter Armor Biliary Split Shield

Biopsy Procedures, Transjugular Intraheptic Portosystemic Shunt (TIPS)



#### Scatter Armor Peripheral Shield

Peripheral procedures where the access is through the femoral Artery.



### Scatter Armor Dialysis Drape


Dialysis Access and Declotting Procedures.











#### **1001 - C Arm Drape**



This item fits all full-size C Arms. Made from durable, high quality polyethylene materials with abrasion resistance to protect from rips or tears.

**HIGHLIGHTS:** CE Certified, Sterile Product, and Latex Free

**QUANTITY:** 30 Drapes (15 drapes per box, 2 boxes included)

 DIMENSIONS:
 42" x 72"

 PACKING SIZE:
 6" x 12" x 17"

 TOTAL WEIGHT:
 13.3 lbs.

#### 1002 - C Arm Kit



Contains C-Arm drape, gray Foot Switch Cover, Image Intensifier and X-ray Tube Cover. Fits OEC® 9600/9800/9900 matches OEC® part #00-901501-01 and E9100AG. Includes 5 mounting clips. Made from durable, high quality polyethylene materials with abrasion resistance to protect from rips or tears.

**HIGHLIGHTS:** Sterile Product, and Latex Free

**QUANTITY:** 20 Drapes per box

**DIMENSIONS:** Multiply Dimensions in this Kit

PACKAGING SIZE: 7" x 21" x 24" TOTAL WEIGHT: 16.7 lbs.

#### 4003 - 54" x 63" Mini C Arm Drape



8003/5406 - Foot Switch Cover



Fits mini-size C Arms, Xi-Scan, Fluoroscan 3" or 4" FOV. Made from durable, high quality polyethylene materials with abrasion resistance to protect from rips or tears.

**HIGHLIGHTS:** Sterile Product and Latex Free

**QUANTITY:** 40 Drapes (20 drapes per box, 2 boxes included)

**DIMENSIONS:** 54" x 63" **PACKAGING SIZE:** 10" x 17" x 24"

TOTAL WEIGHT: 18 lbs.

Foot Switch Cover that is clear and fits with Cinch Lock. Made from durable, high quality polyethylene materials with abrasion resistance to protect from rips or tears.

**HIGHLIGHTS:** Sterile Product, and Latex Free

**QUANTITY:** 100 per box **DIMENSIONS:** 17" x 15" **PACKAGING SIZE:** 10" x 14" x 22"

TOTAL WEIGHT: 8 lbs.


### **Cassette Covers**

### **Protect Your Cassettes**

COMPLETE, HIGHLY EFFECTIVE PROTECTION FROM CONTAMINATION



14" x 17"



Protect your cassettes in our simple and convenient INFAB Cassette Covers.

Available in boxes of 100 or cases of 500 in an easy-to-use dispenser box.



16" x 19"

Cassettes come in contact with blood and other potentially infectious materials.

- Protect your cassettes from cross-contamination
- Protect your patients from infection
- Specially designed to fit your X-ray cassettes
- For use in any imaging area where a risk of contamination exists.

IN-STOCK SAME-DAY SHIPPING



24" x 20"

#### Help Stop the Spread of Infection



The Centers for Disease Control and Prevention (CDC) estimates that roughly 1.7 million hospital-associated infections, from all types of bacteria combined, cause or contribute to 99,000 deaths each year.<sup>1</sup>

<sup>1</sup> Klevens RM, Edwards JR, Richards CL, et al. (2007). "Estimating health care-associated infections and deaths in U.S. hospitals, 2002". Public Health Rep 122 (2): 160–6. PMID 17357358.

<sup>2</sup>Scott II, R. Douglas, (2009). "The Direct Medical Costs of Healthcare-Associated Infections in U.S. Hospitals and the Benefits of Prevention".

National Center for Preparedness, Detection, and Control of Infectious Diseases

#### **Zip Top Cassette Covers**

#### BOX OF 100/500





Zip Top



- Zip top
- Packed in handy dispenser boxes
- 100 or 500 individual covers per box
- · Available in the five cassette sizes
- Always in stock
- Impervious to fluids
- Available in money saving cases of 500.

- · Easy to use
- High quality
- Sanitary
- Disposable
- Top grade polyethylene plastic

#### **FLIP TOP Cassette Covers**

#### **BOX OF 100/500**





Sandwich bag style flap

- "Sandwich" bag style flap
- Packed in handy dispenser boxes
- 100 or 500 individual covers per box
- · Available in the five cassette sizes
- Always in stock
- · Impervious to fluids
- · Available in money saving cases of 500.

- · Easy to use
- High quality
- Sanitary
- Disposable
- Top grade polyethylene plastic

Model	Description		Dimensions	Fits Cassette Size
911417-ZT	Cassette Cover Zip Top	Disposable Box of 100	17" x 19.875"	14" X 17"
901417-ZT	Cassette Cover Zip Top	Disposable Box of 500	17" x 19.875"	14" X 17"
911619-ZT	Cassette Cover Zip Top	Disposable Case of 100	19.25" x 22.25"	16" x 19"
901619-ZT	Cassette Cover Zip Top	Disposable Case of 500	19.25" x 22.25"	16" x 19"
912420-ZT	Cassette Cover Zip Top	Disposable Case of 100	20" x 24.625"	24" x 20"
902420-ZT	Cassette Cover Zip Top	Disposable Case of 500	20" x 24.625"	24" x 20"

Model	Description		Dimensions	Fits Cassette Size
911417	Cassette Cover	Disposable Box of 100	17" x 19.875"	14" X 17"
901417	Cassette Cover	Disposable Case of 500	17" x 19.875"	14" X 17"
911619	Cassette Cover	Disposable Box of 100	19.25" x 22.25"	16" x 19"
901619	Cassette Cover	Disposable Case of 500	19.25" x 22.25"	16" x 19"
912420	Cassette Cover	Disposable Case of 100	20" x 24.625"	24" x 20"
902420	Cassette Cover	Disposable Case of 500	20" x 24.625"	24" x 20"



#### **PRESCRIPTION**

We offer plano, single vision, bifocal, and bifocal progressive options.

#### **PROTECTION**

All glasses have lead protection.
\* Varies per model

#### **FRAMES**

Wide variety of frames from top designers for all face shapes.

#### **PERSONALIZE**

Engrave a name, initials, phone number, or date.



Lens Width: 54mm

Bridge Width: 20mm

Temple Length: 135mm









#### **Oakley Crosslink Zero**

The Oakley Crosslink Zero lead glasses have fixed temples with a three-Point Fit to ensure lenses are aligned for visual correction. Wire cores hidden inside the temples of this lightweight frame offer adjustability for a perfect fit.

- · Rubberized Nosepads and Ear socks for comfort
- Lightweight O Matter™ Frame

FRAME SIZE: Universal WEIGHT: 61.6g

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Black PROTECTION: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case



### Oakley Jupiter Squared

The sporty Jupiter Squared is attitude with refinement. The handles of switchblade knives inspired the earstems and their metal bands and rivets. The Matte Black frames of the Flag Collection are highlighted with American flag emblems.

- · Rubberized Earsocks for comfort and grip
- · USA flag logo on temple arms
- Lightweight O Matter<sup>™</sup> Frame

FRAME SIZE: Universal **WEIGHT:** 74.1q

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Black PROTECTION: .75mm LE

Anti-fog spray, safety strap add-on, protective case INCLUDES:



Lens Width: 56mm

Bridge Width: 18mm

Temple Length: 135mm











Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit. Core wire memory metal temple tips for on the fly eyewear adjustability and grip. Lightweight, durable nylon frame

**FRAME SIZE:** Universal **WEIGHT:** 57G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Matte Dark Grey - Black, Matte Gridiron - Pacific Blue, Matte Black - Space Blue

**PROTECTION:** .75mm LE

INCLUDES: Anti-fog spray, safety strap add-on, protective case

Lens Width: 56.1mm Bridge Width: 15mm Temple Length: 140mm Lens Height: 38.8mm









Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit. Core wire memory metal temple tips for on the fly eyewear adjustability and grip. Lightweight, durable nylon frame

FRAME SIZE: Universal WEIGHT: 60G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive **COLORS:** Matte Black - Gym Red, Matte Midnight Navy - Obsidian

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 54mm Bridge Width: 16mm Temple Length: 140mm Lens Height: 38.9mm















Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit. Core wire memory metal temple tips for on the fly eyewear adjustability and grip. Lightweight, durable nylon frame

FRAME SIZE: Universal WEIGHT: 54G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Dark Grey - Volt, Matte Black, Matte Midnight Navy - Wolf Grey

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 55mm Bridge Width: 14mm Temple Length: 140mm Lens Height: 36.5mm







### ✓ NIKE - 7284

Nike 7284 is a full rim, progressive eligible frame for men and women made of injected plastic. This model features a rounded edges shape, a single bridge, and durable hinges.

FRAME SIZE: Universal WEIGHT: 57G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Brown Basalt - Clear Fade, Clear - Cool Grey

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Bridge Width: 16mm Temple: 145mm Lens Width: 54mm Lens Height: 39mm















Nike 7286 is a full rim, progressive eligible frame for men and women made of eco-friendly plastic, bio acetate. This model features a rectangular shape with a single bridge featuring durable hinges.

**FRAME SIZE:** Universal **WEIGHT:** 68G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Brown Basalt, Matte Obsidian - Blue

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 56.1mm Bridge Width: 15mm

Temple Length: 140mm

Lens Height: 38.8mm









### \_

### NIKE - Adrenaline 22

The Adrenaline 22 Nike frame feature a ventilated rubber nose bridge to reduce fogging. In addition, secure wrap temples provide stability, comfort, and added grip. Durable, no-screw hinges stay securely open and closed.

**FRAME SIZE:** Universal **WEIGHT:** 68G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Black - Volt, Matte Black - Grey

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 64m

Bridge Width: 15mm

Temple Length: 130mm

Lens Height: 40.4mm















The Nike Brazen glasses combine streamlined style, comfort and lightweight durability. They feature a nylon frame and ventilated nose bridge. The Brazen is our most popular model.

FRAME SIZE: Universal **WEIGHT:** 77G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Black- Military Blue, Shiny Black

PROTECTION: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 60.5mm Bridge Width: 16mm Temple Length: 135mm Lens Height: 37.5mm







### ✓ NIKE - Brazen Fury

Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit. Core wire memory metal temple tips for on the fly eyewear adjustability and grip. Lightweight, durable nylon frame

FRAME SIZE: Universal **WEIGHT:** 82G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Dark Grey - Silver Flash, Matte Midnight Navy,

Matte Black - White Red

PROTECTION: .75mm LE

Anti-fog spray, safety strap add-on, protective case **INCLUDES:** 

Lens Width: 61m Bridge Width: 18mm Temple Length: 130mm










Lens Height: 41mm









**BRAZENFURY-03** Matte Midnight Navy





### **NIKE - Brazen Shadow**

The Nike Brazen Shadow features a rubber nose pad and temple tip inlays designed to become tackier as the skin perspires for a solid and secure eyewear fit. Core wire memory metal temple tips provide on-the-fly eyewear adjustability and grip.

FRAME SIZE: Universal WEIGHT: 60G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Anthracite - Grey, Matte Dark Grey - Red

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 59mm Bridge Width: 16mm Temple Length: 130mm Lens Height: 37.4mm









### NIKE - Maverick

The Nike Maverick glasses are pure performance on the inside. The lightweight frame features ventilation and temple arms that provide a soft, firm grip, ensuring you can wear these for hours at a time.

**FRAME SIZE:** Universal **WEIGHT:** 78G

 PRESCRIPTION:
 Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

 COLORS:
 Matte Midnight Navy, Matte Cool Grey, Matte Black Red

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 59m Bridge Width: 15mm Temple Length: 145mm Lens Height: 45.3mm













### **NIKE - Maverick RGE**

The Nike Maverick RGE sunglasses continue the Maverick legacy with long-wearing comfort and a flexible fit in an updated shape. Full rubber nosebridge and temple tips designed to become tackier as the skin perspires for a strong and secure fit.

**FRAME SIZE:** Universal **WEIGHT:** 78G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Matte Midnight Navy- Silver Flash, Matte Black-Grey, Matte Black-Green

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 59mm

Bridge Width: 15mm

Temple Length: 145mm

Lens Height: 40.4mm









### NIKE - Premier

Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure fit.

**FRAME SIZE:** Universal **WEIGHT:** 64G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Matte Black, Matte Anthracite - Grey

Matte Black-Grey

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 60mm

Bridge Width: 13mm

Temple Length: 135mm

Lens Height: 39.3mm

Matte Black-Green

















### NIKE - Rabid 22

Stay clear and focused at every angle in the Nike Rabid 22 lead glasses. Added rubber nose bridge and temples prevent slippage and fogging for optimum performance during procedures.

FRAME SIZE: Universal

**WEIGHT:** 60G

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Matte Black - Red, Matte Obsidian - Blue, Matte Black - White

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 60.5mm Bridge Width: 16mm Temple Length: 135mm







Lens Height: 38mm



### **NIKE - Wave**

The Nike Wave lead glasses feature a slim profile. Textured rubber on temple tips adds stability and grip.

**FRAME SIZE:** Universal **WEIGHT:** 57G

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Ironstone - Grey **PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case







Temple Length: 140mm

Lens Height: 42.6mm











## COMBAT FOG AND GLARE

New specialized lens coating: Anti-Fog or Anti-Reflective

### **DESCRIPTION:**

Lens coatings are available on most glasses and will also increase durability and resistance to scratches.

### **AVAILABILITY**

These coatings are available on available on Plano, Single Vision and Bi-Focal Progressive prescription glasses. \*Only one coating is possible per pair of glasses

**AR** • Anti-Reflective Coating



### Highlights:

- Reduce glare from screens
- Less eye fatigue
- Durable and resistant to scratches

AF · Anti-Fog Coating



### **Highlights:**

- Enhance vision
- Eliminate fog
- · Increase visibility



### GG0752-003

Discover the GG0752-003 Gucci Lead Glasses by INFAB, a harmonious fusion of protection and style. These rectangular Gucci optical glasses feature a transparent grey frame with blue and red stripes, representing a perfect balance of luxury and contemporary fashion.

FRAME SIZE: 56mm - 16mm - 150mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Grey/Blue/Red

WEIGHT: 68 G

**PROTECTION:** .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 56mm

Bridge Width: 16mm

Temple Length: 150mm









# GG0424-001 Black

### GG0424-001

Introducing the GG0424-001 Black Eyeglasses by Gucci, a symbol of style and sophistication. These exquisite eyeglasses boast a classic rectangular full-rim design with a sleek black frame made from a combination of plastic, acetate, and metal. Elevate your look with the timeless elegance and unmistakable luxury of Gucci eyewear.

FRAME SIZE: 56mm - 16mm - 145mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

 COLORS:
 Black

 WEIGHT:
 68 G

 PROTECTION:
 .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 56mm

Bridge Width: 16mm

Temple Length: 145mm







### GG0372-004

Introducing the GG0372-004 Blue/Pink Eyeglasses by Gucci, a captivating blend of style and sophistication. With a mesmerizing blue front and temple color, these eyeglasses exude elegance, while the clear lenses add a fashionable and functional touch.

FRAME SIZE: 51mm - 17mm - 140mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Blue/Pink
WEIGHT: 68 G

PROTECTION: .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 51mm

Bridge Width: 17mm

Temple Length: 140mm









# GG0343-007 Multicolor

### GG0343-007

Sylvie and Web concept glorified by the tri-colored temples with the House tones and new ones. Easy to wear rectangular frame characterized by dark shades in strong contrast with the colored temples. Soft spoken branding with golden Gucci lettering logo on the temple.

**FRAME SIZE:** 57mm - 18mm - 150mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Multicolor
WEIGHT: 79 G
PROTECTION: .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 57mm Bi

Bridge Width: 18mm Temple Length: 150mm

Lens Height: 40mm









### GG0692-004

Modern, light and high-quality surgical stainless steel Gucci frame, rectangular masculine shape, and Gucci logo on temples.

**FRAME SIZE:** 57mm - 16mm - 150mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Black
WEIGHT: 83 G

**PROTECTION:** .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 57mm Bridge Width: 16mm Temple Length: 150mm Lens Height: 38.2mm











### GG0560-005

Vented radiation protection goggle with neoprene foam on the inside designed to protect against fluid splashing and perspiration. Easy to adjust with elastic head band. Prescription compatible (for most Rx). Lightweight and comfortable. Large lens area for maximum protection.

**FRAME SIZE:** 55mm - 20mm - 145mm

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Black/Crystal

WEIGHT: 73 G

PROTECTION: .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 55mm Bridge Width: 20mm Temple Length: 145mm Lens Height: 41.6mm









### **GUCCI GLASSES**

### GG0640-002

Gucci Collection Iconic shapes interpreted in a easy to wear and elegant mood. Rectangular frame. Gucci interlocking logo on temples framed by two gold-coloured bars. Classic color palette.

FRAME SIZE: 55mm - 16mm - 145mm

Alternative Fit

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Havana
WEIGHT: 73 G

**PROTECTION:** .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 55mm Bridge Width: 16mm Temple Length: 145mm Lens Height: 38.7mm











### GG1046-004

Modern, light and high quality Gucci frame with exclusive spring hinge, rectangular shape, and Gucci logo on temples.

**FRAME SIZE:** 55mm - 18mm - 150mm

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Black
WEIGHT: 75 G

**PROTECTION:** .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 55mm Bridge Width: 18mm Temple Length: 150mm Lens Height: 38.2mm







### GG0758-003

Best Seller shape interpreted in a modern, light and high quality version. Rectangular shape. Temples with Gucci logo. Classic color palette with solid and transparent tones.

**FRAME SIZE:** 56mm - 17mm - 150mm

Alternative Fit

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Grey
WEIGHT: 65 G

**PROTECTION:** .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 56mm

Bridge Width: 17mm

Temple Length: 150mm

Lens Height: 38.6mm











### GG0758-001

Best Seller shape interpreted in a modern, light and high quality version. Rectangular shape. Temples with Gucci logo. Classic color palette with solid and transparent tones.

**FRAME SIZE:** 56mm - 17mm - 150mm

Alternative Fit

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Black
WEIGHT: 65 G
PROTECTION: .75 mm LE

INCLUDES: Case, Cloth & Anti-Fog Spray

Lens Width: 56mm

Bridge Width: 17mm

Temple Length: 150mm

Lens Height: 38.6mm











### **Incredibles**

The Incredibles are a stylish wrap around safety frame made of TR-90 high impact resistant plastic. ANSI Z87.1+ approved, lightweight and comfortable. Includes a high quality carrying case. Incredibles glasses are .75mm LE and available from stock in non-prescription.

**FRAME SIZE:** Universal **WEIGHT:** 75G

PRESCRIPTION: Plano, Single Vision (-SV), Bifocal (-BF), Bifocal Progressive (-BF-P)

**COLORS:** Shiny Black, Silver

PROTECTION: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 62m

Bridge Width: 18mm

Temple Length: 145mm

Lens Height: 36mm





















830251 Red

830252 Silver

800070 Blue

### Maxx 30

Prescription compatible (for most Rx), First protective eyewear to receive the American Ophthalmology Associations seal of acceptance. Meets ANSI Z87.1+ standards along with ASTM F803 impact resistance standards. Temple bars and nose bridge have soft rubberized padding for extra comfort, lightweight and comfortable, large lens area for maximum protection Bifocal Progressive Rx frame (requires segment height).

FRAME SIZE: Universal WEIGHT: 63G

PRESCRIPTION: Plano, Single Vision (-SV), Bifocal (-BF), Bifocal Progressive (-BF-P)

**COLORS:** Black, Red, Silver, Blue

PROTECTION: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 53mm

Bridge Width: 17mm

Temple Length: 130mm

Lens Height: 35mm









### **Achilles**

Vented radiation protection goggle with neoprene foam on the inside designed to protect against fluid splashing and perspiration. Easy to adjust with elastic head band. Prescription compatible (for most Rx). Lightweight and comfortable. Large lens area for maximum protection.

**FRAME SIZE:** Universal **WEIGHT:** 69G

PRESCRIPTION: Plano, Single Vision (-SV), Bifocal (-BF), Bifocal Progressive (-BF-P)

**COLORS:** Blue, Plated Silver, Shiny Black

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 53mm

Bridge Width: 17mm

Frame width: 128mm

Lens Height: 35mm













Shiny Black







### Zone

Rubber nose pads and rubberized temple bars. Wrap around design. Lightweight and comfortable. Large lens area for maximum protection.

FRAME SIZE: Universal WEIGHT: 75G

**PRESCRIPTION:** Plano, Single Vision (-SV), Bifocal (-BF), Bifocal Progressive (-BF-P)

**COLORS:** Black, Silver, Red **PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 60mm

Bridge Width: 18mm

Temple Length: 135mm

Lens Height: 36mm









### **GLASSES**

### Evo

The Evo includes permanent side shields with an adjustable nosepiece. These maximum protection glasses come with .75mm protection and a stylish, sleek look. Side shields for maximum protection. Clear inside plastic frame. Adjustable nose piece. EVO is an RX-friendly frame, and it will hold an SV, BF and BF-P.

FRAME SIZE: Universal WEIGHT: 90.4G PRESCRIPTION: Plano

COLORS: Matte Black

**PROTECTION:** FRONT: .75mm LE, SIDE: .50mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 55mm

Bridge Width: 17mm

Temple Length: 140mm

Lens Height: 42mm



















### Zion

The Zion features a large protection area including permanent side shields. These heavy duty glasses are a combination of stylish and protective. Wrap around frame style. Protection on the sides provides a maximum protection of .75mm.

FRAME SIZE: Universal
WEIGHT: 91.4G
PRESCRIPTION: Plano
COLORS: Matte Black

**PROTECTION:** FRONT: .75mm LE, SIDE: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 56mm

Bridge Width: 19mm

Temple Length: 135mm

Lens Height: 42mm









### **Fitovers MX**

The Fitover MX fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitover MX's temples can be easily adjusted for a secure, comfortable fit. The Fitover MX is ideal for fitting over wide prescription eyewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the Fitover MX is a popular fitover model. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

**FRAME SIZE:** Fitover **WEIGHT:** 91G

PRESCRIPTION: Plano

COLORS: Matte Black
PROTECTION .75 mm LE

SIDE SHIELDS 0.50mm (standard) Pb Strips

INCLUDES Anti-fog spray, safety strap add-on, protective case

Lens Width: 69mm

Bridge Width: 17mm

Temple Length: 145mm











Designed to be worn over Rx frames not exceeding 142 mm wide x 40 mm high (5 5/8" x 1 5/8")





Designed to be worn over Rx frames not exceeding 142 mm wide x 50 mm high (5 5/8" x 2")

### **Fitovers XL**

The Fitover XL is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitover XL lead glasses' temples can be benched (adjusted) easily for a secure, comfortable fit. The Fitover XL is ideal for fitting over standard prescription eyewear with square or rounded lenses. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

FRAME SIZE: Fitover
WEIGHT: 108G
PRESCRIPTION: Plano

COLORS: Matte Black
PROTECTION .75 mm LE

SIDE SHIELDS 0.50mm (standard) Pb Strips

**INCLUDES** Anti-fog spray, safety strap add-on, protective case

Lens Width: 64

Bridge Width: 18mm

Temple Length: 150mm







### **Fitovers Razor**

Fits frames not exceeding 132mm width x 38mm height. Wrap around Radiation protection, Hypo-allergenic Crystal Nylon™ Frame. Designed to be worn over small, square, or rectangle prescription eyewear frame shapes.

FRAME SIZE: Universal

WEIGHT: 82G
PRESCRIPTION: Plano

COLORS: Matte Black, Gun Metal, Olive-Charcoal

**PROTECTION** .75 mm LE

INCLUDES Anti-fog spray, safety strap add-on, protective case

Total Width: 132mm Lens Height: 38mm











### **Fitovers Ecofit**

Wrap around Radiation protection with .35mm LE side shield protection. The Hypo-allergenic Crystal Nylon™ Frame is designed to be worn over prescription eyewear

FRAME SIZE: Universal

WEIGHT: 90 G
PRESCRIPTION: Plano

COLORS: Black, Tortoise
PROTECTION .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Total Width: 135mm Lens Height: 41mm





### **Fitovers Lotus**

Fits frames not exceeding 141mm width x 40mm height. Wrap around Radiation protection. Hypo-allergenic Crystal Nylon™ Frame. Designed to be worn over medium, oval, or round prescription eyewear frame shapes.

FRAME SIZE: Universal WEIGHT: 84G

PRESCRIPTION: Plano

COLORS: Matte Black, Smoke Marble, Claret Stripe

**PROTECTION** .75mm LE

**INCLUDES** Anti-fog spray, safety strap add-on, protective case

Total Width: 141mm

Lens Height: 40mm









### **Premium Wraps**

The Premium Wraps feature a large prescription compatible front lens area with maximum radiation protection and side lens for peripheral vision. The frames include spring hinges, formable side temple bars and a fully adjustable nose piece.

FRAME SIZE: Universal
WEIGHT: 97.3G
PRESCRIPTION: Plano
COLORS: Tortoise Shell

Tortoice orien

**PROTECTION:** FRONT: .75mm LE, SIDE: .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 57mm

Bridge Width: 19mm

Temple width: 140mm

Lens Height: 59mm









### **Cyber Flex**

High Quality Safety Glasses. Adjustable temple bar. Large viewing area. Side protection, 0.50mm with an oversized protective lens area.

FRAME SIZE: Universal WEIGHT: 95.4G

**PRESCRIPTION:** Single Vision, Plano

**COLORS:** Clear-Blue **PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: 59mm Bridge Width: 19mm Temple Length: 135-150mm Lens Height: 40mm









### Flip Top

"Flip Top" radiation goggles are an all new product, designed to reduce scatter radiation by as much as 70% (actual reduction is dependent upon imaging technique used) Lightweight, Large protective area: 52mm x 154mm. Flips open when not needed. Adjustable nose bridge. Can be used as a stand alone radiation reduction goggle or fits easily over your existing glasses.

FRAME SIZE: Universal
WEIGHT: 44G
PRESCRIPTION: Plano
COLORS: Amber
PROTECTION: .07mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Lens Width: N/A Total Width: 154mm Temple Length: N/A Lens Height: 52mm







**INFAB Radiation Protection Surgical Gloves** 

### Thinnest Lead-Free Gloves for Maximum Safety

Our Radiation Protection Surgical Gloves set the benchmark in combining safety and dexterity. Our lead-free, powder-free, and latex-free gloves are the thinnest available, offering outstanding radiation protection while maintaining the tactile sensitivity essential for surgical precision.

Thickness	.07" / 0.18 mm
Color	Black
Material	Tungsten
Sizes	6.5, 7.0, 7.5, 8.0, 8.5, 9.0
Compliance	Tested to ASTM F2547-06

	Attenuation Propertie	s
Source Potential	Aluminum Half-Value Level	Measured Attenuation
60 kVp	2.9	56.5%
80 kVp	4.0	47.1%
100 kVp	5.1	41.6%
120 kVp	6.2	37.7%

		1 (	
•	Lead	n-t	ree

- Latex-free
- Powder-free
- · Ultra-thin design
- ASTM tested
- · Multiple sizes available

Size	Model
6.5	797165
7	797170
7.5	797175
8	797180
8.5	797185
9	797190





### XGUARD® Gloves

Superior Radiation Protection with Enhanced Sensitivity

XGUARD® Gloves are engineered for procedures where precision is key. These powder-free, anatomically correct thin gloves are designed to provide superior radiation attenuation while offering superb tactile sensitivity.

- Powder-free
- · Reduced-latex formulation
- · Superior tactile sensitivity
- · Enhanced grip texture
- Multiple sizes available
- · Compliant with standards

Thickness	RR-1 - 0.009" / 0.22 mm
	RR-2 - 0.012" / 0.30 mm
Color	Orange
Material	Natural Rubber Latex and Lead
Sizes	6.5, 7.0, 7.5, 8.0, 8.5, 9.0
Compliance	Meeting EN 455, EN 420, EN 61331 standards, and more

Size	RR-1 Model	RR-2 Model
6.5	RR165-B	RR265-B
7	RR170-B	RR270-B
7.5	RR175-B	RR275-B
8	RR180-B	RR280-B
8.5	RR185-B	RR285-B
9	RR190-B	RR290-B

Box of 5 Individually Wrapped

Attenuation Properties EN 61331-1: 2014	Source Potential	Cuff Attenuation	Palm Attenuation	Finger Attenuation
RR-1 (0.22 mm)	60 kVp	45%	51%	50%
	80 kVp	35%	40%	39%
	120 kVp	29%	33%	32%
RR-2 (0.30 mm)	60 kVp	56%	57%	56%
	80 kVp	44%	45%	57%
	120 kVp	37%	38%	38%



## **Neoprene Gloves**

### Optimal Protection with Tactile Sensitivity

Designed to reduce the amount of scattered radiation exposure to the hands from the primary X-ray beam during radiological procedures. These gloves are Latex-Free, Lead-Free, and of course radiation protected.

	Lead-f	ree
-	Leau-i	166

- Latex-free
- Powder-free
- High sensitivity
- 30% reduction of direct beam radiation transmitted at 100 kVp
- Packed in sterile pouch
- · Shipped from stock

Color	White
Sizes	6.5, 7.0, 7.5, 8.0, 8.5, 9.0

#### **Attenuation Properties**

Source Potential	Cuff Attenuation	Finger Attenuation
60 kVp	38%	49%
80 kVp	30%	40%
100 kVp	24%	33%
120 kVp	19%	27%

Size	Model
6.5	683310
7	683315
7.5	683320
8	683325
8.5	683330
9	683335



### **Revolution Radiation Gloves**

### Superior Safety and Enhanced Touch Precision

INFAB's Exclusive Revolution Radiation
Protection gloves are designed to be the most
comfortable, flexible, lightweight and protective
gloves. Revolution Radiation Protection Gloves
come with tactile grip finger tips. The tactile
grip helps with keeping a strong grasp on
instruments during a long procedure.

- Latex-free
- Powder-free
- Tactile grip
- · Tactile sensitivity
- Sterile
- Maximum tactile sensitivity
- Anatomical fit

Color	Black
Sizes	6.5, 7.0, 7.5, 8.0, 8.5

Size	Model Light	Model Medium	Model Heavy
6.5	700165	735065	775065
7	700170	735070	775070
7.5	700175	735075	775075
8	700180	735080	775080
8.5	700185	735085	775085

Attenuation Properties					
Source Potential	Model Light	Model Medium	Model Heavy		
60 kVp	40%	68%	75%		
80 kVp	35%	80%	69%		
100 kVp	30%	54%	63%		



### **Revolution Maxi-Flex Gloves**

### Supreme Comfort and Dexterity

The perfect way to ensure your hands and lower forearms are fully protected by the most flexible, comfortable radiation protection gloves on the market. Comfortable fit with maximum dexterity. Soft removable, washable and replaceable inner liners reduce odors and allow for washing. Durable outer covers can be replaced. Revolution Maxi-Flex 5 Finger Lead Gloves come exclusively in our Revolution fabric, the softest and most comfortable material available for radiation protection.

Models 683300- Gloves (Pair)

683306- Replacement Liners

683307- Replacement Covers

Slze 15"

Color **Revolution Colors** 

Protection .50mm LE

- · Revolution Fabric in solid colors
- Tough, durable outer coverings
- · Soft, flexible and comfortable
- Removable, washable and replaceable "Cool Wear" moisture-wicking inner liners
- · Replaceable outer covers







































711- White









### **Split Palm Mittens**

Outstanding manual dexterity. Slotted palm allowing the fingers to be extended from the mitten. Complete hand and lower arm protection.

**SIZES:** 15"

**COLORS:** Revolution Colors

**PROTECTION:** .50mm LE

**OPTIONS:** 683700- Open Palm Mittens (Pair)

683701- Open Palm Mittens (Right Only) 683702- Open Palm Mittens (Left Only)

### **Split Palm Mittens with Flap**

Special flap to close the opening in the mitten when finger extension is not desired. Outstanding manual dexterity. Slotted palm allowing the fingers to be extended from the mitten. Complete hand and lower arm protection.

**SIZES:** 15"

**COLORS:** Revolution Colors

**PROTECTION:** .50mm LE

**OPTIONS:** 683710- Open Palm Mittens (Pair)

683711- Open Palm Mittens w/Flap (Right Only) 683712- Open Palm Mittens w/Flap (Left Only)







684000 Hand Shields

### Hand Shields - 684000

Durable, comfortable and fully adjustable shields. Full upper hand protection with total mobility of the hand and fingers. Finger retainers are constructed of durable 1" polypropylene webbing. Hook and loop closure at wrist offers plenty of adjustment.

SIZE: One Size Fits All

**COLORS:** Blue

**PROTECTION:** .50mm LE

- Durable, comfortable and fully adjustable
- Total mobility of the hand and fingers
- Hook and loop closure

### **Maxi-Flex Mittens**

The perfect way to ensure your hands and lower forearms are fully protected by the most flexible, comfortable radiation protection mitten on the market.

**SIZE:** 15"

**COLORS:** Revolution Colors

**PROTECTION:** .50mm LE

**OPTIONS:** 683303- MAXI-FLEX Mittens (Pair)

683308- MAXI-FLEX Mittens Replacement Liners

683309- MAXI-FLEX Mittens Outer Cover

- · Soft, flexible and comfortable
- Removable, washable and replaceable "Cool Wear" moisture-wicking inner liners
- Replaceable outer covers
- Removable, washable and replaceable inner liners



### **Barriers**

INFAB high-quality radiation protection barriers are purpose-built to handle the daily rigors of hospital use. These smooth-rolling and easy-maneuvering barriers offer the highest levels of radiation protection in the industry.

### WINDOW

Your choice of 2.0mm LE lead glass, completely encased with acrylic panels to help prevent damage and breakage, or .50mm LE acrylic.

### **CASTERS**

Hospital-grade swivel locking casters for easy mobility and positive and secure placement.

### **FRAME**

Powder coated steel frame. Precision TIG welded for years of trouble-free service completely enclosing the leaded glass panel for increased durability and safety.

### **PROTECTION**

2.0mm LE lead glass or .50mm LE acrylic window1.5mm LE base panel2 0.75mm LE Curtains (1.5mm LE.)

### **COLORS**

Standard white, black or grey. Custom options available. Upcharge applies.

### **CUSTOMIZATION**

- Color: Frames can be custom powder coated in any pantone color (upcharge applies)
- · Size: INFAB can produce barriers to your specification
- Logo: Add your logo to any barrier
- · Stainless Steel: Optional stainless steel frame





### **EVO-460** Flat Pack Barrier

WINDOW: Acrylic .50mm Pb LE

**WINDOW SIZE:** 24" H x 30" W

**OVERALL SIZE:** 72" H x 34½" W x 25" D

FLOOR CURTAIN: Optional

COLOR: Aluminum Finish

**WHEELS:** 4 swivel casters, 2 with brakes

**SHIPPING:** BOX: 41½" x 51½" x 11½", 135 lbs / 61.24 kg

- Lightweight & Durable Aluminum Frame
- Wheels for Easy Mobility
- Wheel Locking System for Stability
- Compact Shipping Box
- Cost-Efficient Shipping
- Easy 3-Step Assembly

### **STANDARD BARRIERS**

Available in a variety of heights, widths, and window configurations, INFAB barriers are fully customizable, to any specifications. Custom colors, logos and 100% stainless steel options are available. Full head-to-toe protection with our exclusive base curtain included on every barrier. Your choice of 2.0mm LE lead glass or .50mm LE acrylic window. All bases are 1.50mm LE.



683492

### **Large Window Mobile Barrier**

SIZE: 80¼" H x 40½" W X 25" D **WINDOW:** 72" H x 36" W GLASS: Lead or Acrylic (-AC)

**SHIPPING** Crate: 30" x 78" x 48" 480lbs / 217.72 kg

683476

### X-Tall **Mobile Barrier**

SIZE: 78¼" H x 40½" W X 25" D **WINDOW:** 60" H x 36" W GLASS:

Lead or Acrylic (-AC) **SHIPPING:** Crate: 30" x 78" x 48" 450lbs / 204.12 kg

683488

### X-Wide **Mobile Barrier**

SIZE: 74½" H x 76½" W X 25" D WINDOW: 36" H x 72" W GLASS: Lead or Acrylic (-AC)

SHIPPING: Crate: 32" x 85" x 85" 600lbs / 272.16 kg

683473

### **Mid Tall Mobile Barrier**

SIZE: 741/2" H x 341/2" W X 25" D

WINDOW: 60" H x 30" W GLASS: Lead or Acrylic (-AC) **SHIPPING:** Crate: 73" x 30" x 43"

380lbs/172.36ka

### **BARRIER FEATURES**

- All INFAB Barriers include four hospital grade casters, 2 with locking brakes.
- Fully enclosed glass with our exclusive concussion panels to protect the lead glass.
- Easy-clean Sintra plastic base and high quality powder coated steel frames.



683482

### Wide **Mobile Barrier**

SIZE: 72" H x 521/2" W X 25" D **WINDOW:** 36" H x 48" W GLASS: Lead or Acrylic (-AC)

**SHIPPING:** Crate: 30" x 72" x 60" 480 lbs / 217.72 kg

683465

### **Small Window Mobile Barrier**

SIZE: 72¼" H x 52½" W X 25" D **WINDOW:** 12" H x 12" W

GLASS: Lead or Acrylic (-AC) SHIPPING: Crate: 30" x 72" x 60" 380 lbs / 217.72 kg

683468

### **Medium Window Mobile Barrier**

SIZE: 74" H x 52½" W X 25" D 12" H x 24" W WINDOW: GLASS: Lead or Acrylic (-AC) SHIPPING: Crate: 73" x 30" x 57" 410 lbs / 185.97 kg

683460

### **Standard Mobile Barrier**

SIZE: 72" H x 341/2" W X 25" D WINDOW: 24" H x 30" W GLASS: Lead or Acrylic (-AC) SHIPPING: Crate: 30" x 72" x 42"

335 lbs / 151.95 kg

683484

### Regular **Mobile Barrier**

SIZE: 641/2" H x 401/2" W X 25" D **WINDOW:** 24" H x 36" W GLASS: Lead or Acrylic (-AC) SHIPPING: Crate: 30" x 65" x 48"

430 lbs / 195.04 kg

683458

### "Shorty" **Mobile Barrier**

SIZE: 45" H x 341/2" W X 25" D WINDOW: 24" H x 30" W GLASS: Lead or Acrylic (-AC) SHIPPING: Crate: 45" x 30" x 42" 220 lbs / 99.79 kg



### 683461 Tilted Mobile Barrier

WINDOW: 2.0mm LE Lead Glass or .50mm LE Acrylic

**WINDOW SIZE:** 24" H x 30" W

**OVERALL SIZE** 70" H x 34¼" W X 25" D

**COLORS:** Standard white, grey or black. Custom options

available.

**FRAME:** Metal, Stainless Steel (Upgrade)

**SHIPPING:** Crate: 73" x 30" x 43", 355 lbs / 161.03 kg

Powder coated steel frame

- Acrylic impact panels, enclosed in a precision TIG welded frame for years of trouble-free service, completely enclosing the leaded glass panel for increased durability and safety
- 2 high-quality hospital grade locking casters
- 1.50mm LE opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Ergonomical viewing area
- Shipping worldwide by ground, air or sea freight

Tilted design for better view





### 076993

### **Technician Protection Mobile Screen**

The Technician Protection Mobile Screen can be easily raised and lowered by hand and stays in position thanks to the perfectly counter-balanced internal springs. The high quality steel construction and stainless springs provide a lifetime of trouble-free use.

**SIZE RAISED:** 73" H x 31" W x 20 ½" D

**SIZE LOWERED:** 44" H x 31" W x 20 ½" D

**WINDOW SIZE:** 37.4" H x 27" W

**COLORS:** Standard white, grey or black. Custom options

available.

**FRAME** Metal, Stainless Steel (Upgrade)

**SHIPPING:** Crate: 31" x 49" x 41", 210 lbs/95.25kg

Retractable Window

- Multiple heights
- · Perfectly counter balanced
- Easily raised and lowered
- 1.0mm LE base
- .50mm LE lead acrylic window
- 2 high-quality hospital grade locking casters
- · Complete with bar to hold controls



New steel frame adds

strength and prevents bending

or warping of viewing window



### **REV-MHAS**

### **Revolution Mobile Hanging Acrylic Shield**

Fully adjustable mobile shield that can be altered to user's preference. This product provides the benefits of an overhead suspension barrier with the mobility to fit into exam rooms of any size. The acrylic barrier can be raised, lowered, and rotated 360° in both vertical and horizontal planes to fit snug against the torso of the patient.

#### **KEY FEATURES:**

- · Mobile unit with superior stability.
- · Gives user full control of shield placement.
- · Pivot ball allows for 360° motion.
- · Durable 0.50mm acrylic shield.
- Small footprint, does not interfere with physician's position at the table.

DIMENSIONS: 79" H x 37" W x 25" D

WINDOW SIZE: 24" H x 30" W

PROTECTION: Lead acrylic window .50mm LE

COLORS: Standard off-white powder coated.
 Other colors available.

FRAME Powder Coated Steel

WHEELS: 4 swivel hospital grade casters (2 with brakes)

Full 360° Pivot ball for precise positioning

Front profile



### **REV-PMC**

#### **Revolution Mobile Curtain**

A unique new design for operating physicians that allows unparalleled protection for close-to-table studies. Sturdy, yet thoughtfully engineered with a small footprint so it can easily be positioned exactly where you want it for maximum scatter radiation protection.

**SIZE:** Adjustable height 43½" - 66"

Base: 34" W x 25" D

**CURTAIN:** 38"H x 30"W .50mm LE lead curtain with grey Revolution fabric

FRAME COLORS: Standard white, grey or black. Custom options available.

CURTAIN COLORS: Standard Revolution grey. Custom options available.

FRAME: Powder Coated Steel, Stainless Steel (Upgrade)

**SHIPPING:** Box: 13" x 42" x 9", 45 lbs / 20.41 kg

Box: 26" x 30" x 8", 30 lbs / 13.60 kg



- Special 4 caster T-bar base design to stabilize the unit
- Smooth easy maneuvering
- Special spring mechanism allows for easy raising, lowering, and automatically holding shield
- .50mm LE double protective curtain
- 2 high-quality hospital grade locking casters



### **REV-SPMC**

### **Revolution Solid Panel Mobile Curtain**

**SIZE:** 56" H x 30" W **CURTAIN:** 30" x 30"

**PROTECTION: REV-SPMC10** 1.0mm LE

REV-SPMC50 .50mm LE

**FRAME COLORS:** Standard white, grey or black. Custom options available. **CURTAIN COLORS:** Standard Revolution grey. Custom options available.

**FRAME:** Powder Coated Steel, Stainless Steel (Upgrade)

**SHIPPING:** Box: 21" x 21" x 10", 18 lbs / 8.17 kg

Box: 16" x 33" x 7", 28 lbs / 12.7 kg


• Heavy-duty hospital-grade purpose-built unit.

- · Adjustable height between 37" and 56"
- 30" wide by 30" tall lead vinyl curtain
- 2 high-quality hospital grade locking casters
- Smooth and easy maneuvering



56"











## **REV-MC**

## **Economy Mobile Curtain**

SIZE: 66" H x 28" W

30" x 30".50mm LE curtain with grey Revolution fabric **CURTAIN:** Standard white, grey or black. Custom options available. FRAME COLORS: Standard Revolution grey. Custom options available. **CURTAIN COLORS:** Powder Coated Steel, Stainless Steel (Upgrade) FRAME:

**SHIPPING:** Box: 24" x 30" x 12", 20 lbs /9.07 kg

Box: 6" x 6" x 6", 27 lbs / 12.25 kg

· Raise and lower shield

· Smooth easy maneuvering

• .50mm LE double protective curtain

• 2 high-quality hospital grade locking casters











**REV-LOCKER-20** 

## Revolution Mobile Apron Lockers

Features digital passcode locks with keypad Set up to 3 administrative and user pass Battery operated lock with backup battery and emergency power supply





Passcode Lock with Keypad







## **REV-LOCKER-5**

## Revolution 5 Hanger Mobile Apron Locker

**ASSEMBLED SIZE:** 65" H x 16" W x 26.5" D

**COLORS:** White, Grey, or Black

Custom colors available.

Ships in standard white unless specified.

**SHIPPING:** Crate: 29" x 20" x 70", 170 lbs / 77.11 kg

Includes 5 hangers Secure up to 5 aprons

## **REV-LOCKER-10**

## Revolution 10 Hanger Mobile Apron Locker

**ASSEMBLED SIZE:** 65" H x 32" W X 26.5" D

**COLORS:** White, Grey, or Black

Custom colors available.

Ships in standard white unless specified.

**SHIPPING:** Crate: 29" x 36" x 70", 205 lbs / 93 kg

Includes 10 hangers Secure up to 10 aprons

## **REV-LOCKER-20**

## Revolution 20 Hanger Mobile Apron Locker

**ASSEMBLED SIZE:** 65" H x 48" W X 26.5" D

**COLORS:** White, Grey, or Black

Custom colors available.

Ships in standard white unless specified.

**SHIPPING:** Pallet: 29" x 52" x 70", 230 lbs / 148.33 kg

Includes 20 hangers Secure up to 20 aprons

# **Apron Racks**

INFAB Apron Racks are purpose-built pieces of hospital equipment, designed to last as long as the facility they are being used in.

## **OPTIONS**

- Mobile
- Wall-Mounted
- Single Apron
- Multi-Apron
- Deluxe
- Economy
- Custom Colors

## HOSPITAL GRADE MATERIALS

- Cold-rolled metal or optional stainless steel frames
- Hospital grade casters, two with locking brakes, for easy mobility and positive placement when parked.

## THOUGHTFUL DESIGN

- Easy to use
- Space-efficient
- · Plentiful storage
- Easy on/off apron convenience









#### 10 Hanger Mobile Apron Valet

**ASSEMBLED SIZE:** 58" H x 26½" W X 25" D

**COLORS:** White/Chrome, Custom Options

**SHIPPING:** Box 1: 55" x 7" x 7", 25 lbs / 11.34 kg

Box 2: 16" x 32" x 14", 34 lbs / 15.42 kg

- Space saving; providing major storage while taking up minimum space.
- · Steel construction with baked enamel finish.
- · 2 high-quality hospital grade locking casters
- 10 Deluxe Custom Apron Hangers (683450) included





#### 20 Hanger Mobile Apron Valet

**ASSEMBLED SIZE:** 58" H x 481/2" W X 25" D

**COLORS:** White/Chrome, Custom Options

**SHIPPING:** Box 1: 55" x 7" x 7", 32 lbs / 14.51 kg

Box 2: 16" x 32" x 14", 43 lbs / 19.50 kg

- Space saving; providing major storage while taking up minimum space.
- · Steel construction with baked enamel finish.
- · 2 high-quality hospital grade locking casters
- 20 Deluxe Custom Apron Hangers included
- Optional: full circle loop hangers (-H)



## 683419

## Mini Mobile Apron Rack with Hangers

**ASSEMBLED SIZE:** 78" H x 27" W X 25" D

**COLORS:** White, Custom Options

**SHIPPING:** Box: 29" x 49" x 8", 59 lbs / 26.76 kg

- · Six apron hangers included
- Heavy gauge cold rolled steel
- Powder coated
- · Attractive, durable, long lasting finish
- · 2 high-quality hospital grade locking casters
- Holds up to six protection aprons, a pair of lead gloves and two pairs of radiation protection glasses







## **Deluxe Mobile Apron Rack**

**ASSEMBLED SIZE:** 53" H x 39" W X 26" D **COLORS:** White, Custom Options

**SHIPPING:** Box 1: 16" x 32" x 14", 25 lbs / 11.34 kg

Box 2: 40" x 29" x 12", 55 lbs / 24.95 kg Box 3: 58" x 17" x 9", 52 lbs / 23.59 kg

- · Super heavy-duty
- 2 high-quality hospital grade locking casters
- Stores up to ten protection aprons
- Optional vest and skirt adapter (Part Number 683000)
- · Holds up to 200 lbs. of aprons

## 683426

## 5 Arm Budget Saver Mobile Apron Rack

ASSEMBLED SIZE: 58" H x 26" W X 30" D

COLORS: White, Custom Options

**SHIPPING:** Box 1: 51" x 7" x 7", 31 lbs / 14.06 kg

Box 2: 16" x 32" x 14", 34 lbs / 15.42 kg Box 3: 16" x 41" x 7", 21bs / 9.53 kg Box 4: 26" x 21" x 2", 12 lbs / 5.44 kg

- Economically priced
- · Heavy gauge steel
- Powder coated
- 2 high-quality hospital grade locking casters
- Stores up to 5 protection aprons
- · Convenient additional storage on base
- Optional vest and skirt adapter (Part Number 683000)

#### 683427

### 10 Arm Budget Saver Mobile Apron Rack


ASSEMBLED SIZE: 58" H x 44½" W x 30" D


COLORS: White, Custom Options

**SHIPPING:** Box 1: 51" x 7" x 7", 36 lbs / 16.33 kg

Box 2: 16" x 32" x 13", 34 lbs / 15.42 kg Box 3: 58" x 16" x 9", 40 lbs / 18.14 kg Box 4: 22" x 44" x 2", 21 lbs / 9.53 kg

- · Economically priced
- Heavy gauge steel
- Powder coated
- 2 high-quality hospital grade locking casters
- Stores up to 10 protection aprons
- · Convenient additional storage on base
- Optional vest and skirt adapter (Part Number 683000)





5 Arm Left Folding Apron Rack

683100

5 Arm Right folding Apron Rack

#### 683000

Vest and Skirt Rack Adapter

• Skirt rack adapter easily mounts to rack with included tool. Skirt hanging adapter provides convenient storage for vest & skirt aprons.

#### Wall-Mounted 5 Arm Apron Rack

**ASSEMBLED SIZE:** 8¼ H x 6" W X 23" D **COLORS:** White, Custom Options

- Stores up to 5 aprons
- Swings out easily for immediate access to aprons
- Heavy duty
- Power coated steel mounting bracket and storage arms
- Arms with nylon bushings
- · Racks fold flush to wall
- Can be mounted side by side
- Optional vest and skirt adapter (Part Number 683000)

#### **SHIPPING**

 683100 5 Arm Folding Right Hand Folding
 Box 1: 10 lbs (12" x 30" x 6")

 683101 5 Arm Folding Left Hand Folding
 Box 1: 10 lbs (12" x 30" x 6")



#### **Wall-Mounted Wall Tower Rack**

**ASSEMBLED SIZE:**  $48" \text{ H x } 13 \frac{3}{4}" \text{ W X } 27" \text{ D (w/racks)}$ 

**COLORS:** White, Custom Options

• Includes right and left 5 arm apron racks


- Heavy duty
- Handles toughest hospital environments
- Heavy gauge steel
- Powder coated
- Scratch resistant
- Easy to install
- · Hardware provided
- Requires only minimal space
- 13.75" wide x 48" tall

#### **SHIPPING**

683115 Weak Wall Support Box 1: 10 lbs (12" x 30" x 6") 683415 Rack/Tower W/ R.H. & L.H. Racks Box 1: 10 lbs (12" x 30" x 6")



**683115** Weak Wall Support



#### **Wall Mounted Rack**

**ASSEMBLED SIZE:** 12" H x 36" W X 26" D **COLORS:** White, Custom Options

- · Heavy gauge steel
- 200 lbs. capacity
- Powder coated
- · Scratch resistant
- · Easy to install
- · Requires only minimal space
- Stores up to 10 aprons
- Optional vest and skirt adapter (Part Number 683000)

#### 683429 SHIPPING

Box 1: 40 lbs (16" x 58" x 9")

## 683430

## **Apron and Glove Combination Rack**


ASSEMBLED SIZE: 17" H x 2834" W (91/2" Baseplate)

**COLORS:** White, Custom Options

- Rounded corners to prevent accidental apron tears or skin punctures
- Quick, convenient storage
- Stores 1 apron and 1 pair of lead gloves
- · Sturdy & compact
- Mounts quickly and easily
- · Heavy gauge steel mounting plate
- · Half inch outside diameter tubing
- · Electro static powder coated

#### 683430 SHIPPING

Box 1: 9 lbs. (19" x 30" x 11")



#### 683420

## **Apron and Glove Storage Rack**

**ASSEMBLED SIZE:** 25" H x 28" W x 6" D

**COLORS:** White, Custom Options

- · Quick, convenient storage
- Sturdy
- Compact
- Mounts quickly and easily
- Heavy gauge steel mounting plate
- Half inch outside diameter tubing
- · Electro static powder coated
- Great for dental aprons

#### 683420 SHIPPING

Box 1: 15 lbs (26" x 31" x 8")



#### **Protective Accessories Rack**



ASSEMBLED SIZE: 21" H x 18½" x 5" D

COLORS: White, Custom Options

- Storage unit for two or more aprons, a pair of glasses and a pair of lead gloves
- Mounts quickly and easily
- 18" wide
- Extends a mere 4.5" from the wall
- Contains no sharp edges
- Perfectly suited to small areas
- Powder coated heavy steel construction

#### 683102 SHIPPING

Box 1: 8 lbs (19" x 23" x 8")

## 683424 | Single Bar Apron Rack



**ASSEMBLED SIZE:** 6" H x 18 ½" W x 3" D

COLORS: White, Custom Options

- · Storage unit for one apron
- · Mounts quickly and easily
- Contains no sharp edges
- Perfectly suited to small areas
- Powder coated heavy steel construction

#### 683102 SHIPPING

Box 1: 3 lbs (24" x 4" x 8")

## 683452 | Locker Size Hanger



SIZE:

6" H x 14 ½" W

**COLORS:** 

White, Custom Options



Optional hook (-H)

- Securely hold a radiation protection apron
- Constructed from heavy gauge steel tube
- · Designed to fit in personal locker
- Powder coated

## 683450

## **Deluxe Custom Hanger**



SIZE:

7½" x 22½"

**COLORS:** 

Chrome, Custom Options



Optional hook (-H)

- Securely hold a radiation protection apron
- · Constructed from heavy gauge steel tube
- · Chrome plated
- Available in a pack of 5 hangers
- Part Number 683451

## 683412 Two Peg Wall Mounted Rack



683412 683412-LS

**ASSEMBLED SIZE:** 9" H x 15 %" W x 7%" D **COLORS:** White, Custom Options

**SHIPPING:** Box: 11" x 18" x 9", 6 lbs / 2.72 kg

683411 Eight Peg Wall Mounted Rack



683411

683411-LS

Eight Peg Wall Mounted Rack

Eight Peg Wall Mounted Rack with Locking System

**ASSEMBLED SIZE:**  $9" \text{ H x } 46 \frac{1}{2}" \text{ W x } 7\frac{1}{4}" \text{ D}$  **COLORS:** White, Custom Options

**SHIPPING:** Box: 13" x 47" x 10", 21 lbs / 9.53 kg

683410

**Six Peg Wall Mounted Rack** 



683410 683410-LS

Six Peg Wall Mounted Rack Six Peg Wall Mounted Rack with Locking System

**ASSEMBLED SIZE:**  $9" H x 34 \frac{1}{2}" W x 7 \frac{1}{4}" D$  **COLORS:** White, Custom Options

**SHIPPING:** Box: 11" x 37" x 9", 16 lbs / 7.26 kg

LS Locking System



- Available in 2 peg, 6 peg, or 8 peg configurations
- · Inexpensive and very effective
- · Heavy gauge steel plate
- · Electrostatic powder coated
- Tilted storage bars so aprons wont slip off
- · Requires minimal space
- · Includes lock



## **Curtains**

INFAB offers customizable radiation protection curtains of superior quality, available in various protection levels, sizes, colors, and fabrics. Our durable curtain tracks, tailored to your specifications, come in both straight and curved options for the perfect room fit.

Contact us with your specifications for a customized quote.

#### **Lead Curtain**



Our lead curtains include one grommet hole with a "S" hook every six inches.

COLORS: Classic Nylon Colors

**CORE:** Lightweight Lead, KIARMOR Lead-Free

**PROTECTION:** .25mm LE, .35mm LE, .50mm LE,

.75mm LE, 1.0mm LE

**SIZE:** Per square foot

**683650** .25mm LE curtain with grommets (per sq. ft.)

**683655** .35mm LE curtain with grommets (per sq. ft.)

**683660** .50mm LE curtain with grommets (per sq. ft.)

683670 .75mm LE curtain with grommets (per sq. ft.)

683675 1.0mm LE curtain with grommets (per sq. ft.)

## **Straight Track for Curtain**

Mill finish 6063-T5 extruded aluminum 1-3/4" wide x 2" high, .080 thick. **Priced per linear foot.** 

**INCLUDES:** 1 Straight Track

2 End Stops

1 Single Carrier for each grommet hole on curtain

**OPTIONS:** Clamp Hanger

Ceiling Splice Trim Chain

Single Carrier

2 polyethylene wheels riveted to plated steel body with swivel hook. Straight Track



#### **End Stop**

Bolt with horizontal pipe spacer to prevent carriers from traveling beyond desired location. For use with walk-along tracks.


#### **OPTIONS:**



#### Clamp Hanger

Formed steel brackets clamp to track. Adjusts to any position.

Suspend track with welded chain or wire rope. 1 pair pictured.



#### Ceiling Splice

To align and lock track sections together on ceiling mounted tracks.

1 pair pictured.

Note: Cannot be used with No. 132 Overlap Master or No. 10 Pivot Device due to interference with the wheels.



#### Trim Chain

3 links No.10 single jack chain added to carriers to provide 3" of curtain height adjustment





End Stop

Bolt with horizontal pipe spacer to prevent carriers from traveling beyond desired location. For use with walk-along tracks.

## 683690

## **Curved Track for Curtain**

Mill finish 6063-T5 extruded aluminum 5/8" wide x 1-11/16" high, 0094" thick. For curved track applications only. **Each piece requires custom quoting.** 

**INCLUDES:** 1 Curved Track

2 End Stops

1 Single Carrier for each grommet hole of curtain

**OPTIONS:** Clamp Hanger

Ceiling Splice
Ceiling Hanger
Trim Chain

#### Single Carrier

2 polyethylene wheels riveted to plated steel body with swivel hook.




#### Trim Chain

3 links No.10 single jack chain added to carriers to provide 3" of curtain height adjustment

#### **OPTIONS:**




Splice
Steel plates align and lock track
together. 1 pair pictured.



#### Clamp Hanger

Formed steel brackets clamp to track. Adjusts to any position. Suspend track with welded chain or wire rope.1 pair pictured.



Ceiling Hanger

Steel plate mounts flush to ceiling. Track clamps to bracket with offset clips.

# **Lead Curtains**

These scatter shields are attached to the table to protect both the operator and the patient from harmful radiation.

# **Easy Setup & Installation**

Easily slides onto most tables with side rails and quickly locks in to place with an adjustable knob.

#### 077036

#### **Single Panel Scatter Shield**

The below-table protection screen is attached to the table to protect the operator's lower body from radiation. The upper part of the screen (removable) protects against scatter radiation from the patient's body. The lower part protects the operator's lower body against the radiation from the tube source.


#### **UPPER SCREEN**


- .50mm LE with stiff support
- •11" H x 26" W
- · Hooks into mobile arm that can be easily removed

#### **LOWER SCREEN**

- .50mm LE made up of 4 strips
- 30" L x 26" W

Universal adapter bracket fits almost every table with a side rail





## **Multiple Panel Scatter Shield**

The below-table protection screens are attached to the table to protect the operator's lower body from radiation. The upper part of the screen (removable) protects against scatter radiation from the patient's body. The lower part protects the operator's lower body against the radiation from the tube source.

#### **UPPER SCREEN**

- .50mm LE with stiff support
- 11" H x 26" W
- · Hooks into mobile arm that can be easily removed

#### **LOWER SCREEN**

- .50mm LE made up of 4 strips
- 30" L x 26" W

#### **ADDITIONAL STRIP (INCLUDED)**

- Applied with addition of intermediate joint
- 30" L x 10" W






Available


#### **Saddle Shield Radiation Protection Curtains**

Saddle Shield Radiation Protection Curtians designed to help protect the body from being exposed to radiation.

- .50mm LE
- Significantly reduces scatter radiation
- Overlapping split front panels
- No special hardware or modifications to imaging equipment required
- · Available in any of our fabric colors
- Provide exact dimensions and drawing to fit your specific size table







#### **Patient Protection Set**

**INCLUDES:** Storage rack with 4 aprons (Toddler, Child, Average, Large)

COLORS: Classic Nylon Colors

CORE: Lightweight Lead, KIARMOR Lead-Free
PROTECTION: .25mm LE, .35mm LE, .50mm LE

SIZES: Toddler, Child, Average Adult, Large Adult









683401 683402 9" H x 11.5" W 13" H x 13" W Toddler Apron A

683403 14" H x 17.5" W Average Adult Apron

683404 17" H x 21" W Large Adult Apron



### **Half Apron**

#### **HSV**

Hook and Loop Closure Half Apron

**HSV-S1** 12" x 24" Half apron with hook and loop closure

**HSV-M1** 18" x 24" Half apron with hook and loop closure

**HSV-L1** 24" x 24" Half apron with hook and loop closure



Standard hook and loop closure

#### **ASB**

Buckle Closure Half Apron

**ASB-S1** 12" x 24" Half apron with buckle closure

**ASB-M1** 18" x 24" Half apron with buckle closure

**ASB-L1** 24" x 24" Half apron with buckle closure



#### **Disposable Eye Shields** PP-ES



- Minimum 40% dose reduction
- Box of 50
- · Barium sulfate shields are non-latex and 100% disposable



## **Disposable Thyroid Shield**



- Minimum 60% dose reduction
- Box of 50
- · Barium sulfate shields are non-latex and 100% disposable

#### 684120

#### **Ovarian Shields**



- - Set of the three (Infant, Youth and Adult)
  - Designed especially for ovarian procedures
  - 1.0mm LE
  - Easy to use and explain to your patient
  - · Available in our Classic Nylon Colors
  - Durable
  - · Easily wipes clean

## **Gonad Shields**



- - Set of three or available individually: (Infant, Youth and Adult)
  - 3/4" tie strap
  - Easy fastening around the waist
  - .50mm LE
  - · Available in our Classic Nylon Colors



## **Triangular Gonad Shields**



- Set of three (Infant, Youth and Adult)
- 1/2" inch tie strap
- Easy fastening around the waist
- .50 mm Pb equivalent



## PTD | Patient Transfer Device

- Radiolucent plastic
- 18" or 22" wide models
- Also available with storage rack
- Large hand holds
- One inch wide wrist straps
- Anti-static
- Thin
- Board slides easily from surface to surface
- Tested up to 300 lbs.


MODELS: PTD-18, PTD-22

COLORS: White SHIPPING: PTD-18

Box: 73" x 23" x 1", 11 lbs / 4.99 kg

PTD-22

Box: 73" x 23" x 1", 13 lbs / 5.90 kg



## 683399 Patient Transfer Rack

Wall mounted rack for PTD storage

**SIZE** 9" H X 15" ¼ W x 6 ¾" D

## **Coated Sponges**

Coated Sponges combine medical-grade durability with fluid resistance, ensuring easy cleaning and maintenance. These sponges, designed with a patented flared angle, reduce x-ray imaging artifacts, enhancing image clarity without affecting patient positioning.

- · Durable medical-grade coating
- · Reduced imaging artifacts
- · Easy-to-clean, fluid-resistant
- · Designed for clear, comfortable positioning



**INF-NCAA** 

**Adult Oblique Finger Block - Coated** 

Dimensions 10.5" x 6" x 3"



**INF-NCAD** 

Adult Lateral Finger Block - Coated


Dimensions 13" x 8" x 3"



**INF-NCAK** 

**Circular Head Block - Coated** 

Dimensions Diameter: 7" (top) - 9" (base) x 1.5"



INF-NSGP

## Positioning Sponge General Kit P - Coated

#### Includes:

- INF-NCAA (x1) Adult Oblique Finger Block 10.5" x 6" x 3"
- INF-NCAK (x1) Circular Head Block 7" with 9" base
- INF-NCBD (x1) 45° Spinal Wedge 27.75" x 10" x 7"
- INF-NCBI (x2) 15° Wedge 10" x 12.5" x 3"

INF-NCBD

45 Degree Spinal Wedge - Coated

Dimensions (LxWxH) 27.75" x 10" x 7"





**INF-NCFD** 

45 Degree Spinal Body Wedge Sponge - Coated

Dimensions (LxWxH) 21" x 7" x 7"

## **Closed-Cell Sponges**

Closed Cell Sponges are made from a lightweight, radiolucent, non-porous foam. These positioning sponges are entirely fluid resistant, easy to clean and extremely durable. Featured in a noncoated, premium charcoal grey color.

- · 100% Radiolucent
- · Artifact free imaging
- · Cleanable, fluid resistant surface
- · Lightweight and comfortable



**INF-CCAA** 

**Finger Block Positioner - Closed Cell** 

Dimensions (LxWxH) 5" x 5" x 3"



**INF-CCBD** 

45 Degree Spinal Body Wedge - Closed Cell

Dimensions (LxWxH) 21" x 7" x 7"



**INF-CCBL** 

**Rectangle Block Positioner - Closed Cell** 

Dimensions (LxWxH) 10" x 7.5" x 3"



INF-CCK1

General Sponge Positioning Kit 1 - Closed Cell

#### Includes:

- INF-CCBD (1) 45° Spinal Body Wedge 21" x 7" x 7"
- INF-CCAA (1) Finger Block Positioner 5" x 5" x 3"
- INF-CCAK (1) Circular Head Block 7" x 1.5"
- INF-CCBI (2) 15° Wedge 11" x 7" x 3"
- INF-CCBV (1) Square Block 6" x 6" x 2"
- INF-CCBL (1) Rectangle Block 10" x 7.5" x 3"
- INF-CCMB (1) Rectangle Block 10" x 8" x 2"



INF-CCBN

**Closed Cell Rectangle** 

Dimensions (LxWxH) 18" x 14" x 3"

## **Weight-Bearing Studies**





INF-WBCH

Weight Bearing Cassette / DR Panel Holder

This holder, crafted from high-density crosslink foam, provides a sturdy and cleanable surface for foot and ankle imaging, supporting up to 300 lbs.

Dimensions Foam Size (LxWxH) 17.5" x 11.5" x 4"

(LxWxH) Slit Size (LxWxH) 11.5" x 0.625" x 3"





INF-WBCH2

2 Step Weight Bearing Cassette DR Panel Holder

This two-tiered step holder, made from durable, antimicrobial crosslink foam, supports up to 300 lbs for stable, handrail-free weight-bearing exams with easy maintenance.

Dimensions (LxWxH)

Foam Size: 24" x 16" x 12"

First Step: 10" x 16" x 6" Top Platform: 14" x 16" x 12"

Dimension of the slot: 14"L x 0.625"W x 1.5"D

## **Table Pads**



INF-PPAD-3080

Premium Vinyl Table Pad 30x80x2

This premium table pad offers a 2-inch thick cushion for patient comfort, covered in fire-retardant, antimicrobial vinyl that's easy to clean and ideal for x-ray procedures.

Color

Blue

Dimensions

30" x 80" x 2"



INF-SPAD-2472

Standard Table Pad - 24x72x1

Color

Blue

**Dimensions** 

24" x 72" x 1"

INF-SPAD-3080

Standard Table Pad - 30x80x1

Color

Blue

Dimensions

30" x 80" x 1"

These standard table pads feature a fire-retardant, antimicrobial vinyl cover, providing a secure and comfortable patient surface. They're easy to clean, radiolucent, and designed for durability, fitting a range of imaging equipment.

## **FDA**

Establishment Registration & Device Listing: INFAB, LLC.
1040 Avenida Acaso
Camarillo, CA 93012
Registration Number: 2024351

Date of Registration Status: 2017 Owner/Operator: INFAB, LLC. 1040 Avenida Acaso Camarillo, CA 93012 Owner/Operator Number: 2024351 1040 Avenida Acaso Camarillo, CA 93012 Phone: 805-987-5255

#### **FDA Certificate**

Device Classification Name:	Shield, protective, personnel
510(k) Number	K896289
Device Name	RADIATION PROTECTION APRON
Applicant	INFAB, LLC. 1040 Avenida Acaso Camarillo, CA 93012
Regulation Number	892.6500
Classification Product Code	KPY
Date Received	10/31/1989
Decision Date	12/14/1989
Decision	Substantially equivalent (SE)
Classification Advisory Committee	Radiology
Review Advisory Committee	Radiology
Statement/Summary/Purged Status	Purged, no summary or statement
Туре	Traditional
Reviewed by Third Party	No

## 510(k) Premarket Notification

Proprietary Name	Non leaded apron
Classification Name	Apron, non leaded
Product Code	EAJ
Device Class	1
Regulation Number	892.6500
Medical Specialty	Radiology
Registered Establishment Name	INFAB, LLC.
Registered Establishment Number	2024351
Owner/Operator	INFAB, LLC.
Owner/Operator Number	2024351
Review Advisory Committee	Radiology
Establishment Operations:	Manufacturer













1. PRICING, DELIVERY AND PAYMENT. The prices for Product delivered under this Order do not include customs charges or duties, or sales, use, excise or similar taxes, or special packaging, or transportation, rigging, dravage, insurance charges and other costs associated with the delivery of Product to the destinations specified by the Customer unless specifically set out on the invoice as individual line items. Therefore Customer will pay directly to Seller ("Seller"), in addition to the prices specified in this Order, any amount set forth on Seller's invoices for customs charges, duties and/or taxes (or provide applicable tax exemption certificates acceptable to the taxing authorities) and special packaging costs, and shall pay directly to the involved carrier(s) all such delivery related costs. All licenses or certificates required for import or export of the Product shall be the responsibility of and obtained at the cost of Customer. Seller's terms are cash upon delivery, unless otherwise set forth in a distribution or reseller agreement. If payment is by credit card, an additional fee will apply. If credit terms are approved by Seller. Customer shall pay all invoices net thirty (30) days after date of shipment. If payment is not made when due, Customer shall be assessed a late charge equal to one and one-half percent (1-1/2%) of the past due balance for each month or partial month following the due date. Additionally, at sellers discretion, all discounts shall be removed and any amounts due shall be the current list price of such items, regardless of any other agreement or discount offered. All deliveries shall be made F.O.B. common carrier. Seller's plant. Title and risk of loss or damage will transfer from Seller to Customer upon delivery of goods at the F.O.B. point. Customer shall notify Seller, in writing, within one (1) calendar day after delivery, of any claimed shortages, damage or of rejection of the delivery. Sufficient detail to establish grounds shall be included in any rejection.

2. PRODUCT SAMPLES. Product samples will be provided to licensed distributors of Seller's product. Upon termination of a distribution arrangement, distributor shall return to Seller, at distributors expense, all samples previously provided to distributor by Seller in resalable as new condition. For any goods not returned in accordance with the above, Distributor will be invoiced for sample product at the same discount established for distributor for production goods. If invoice is not paid within 10 days of receipt, all discounts will be canceled and prices for samples will revert to full List Price.

3. SECURITY INTEREST. Seller reserves, until full payment is received, a Purchase Money Security Interest in each Product Unit and the proceeds of any sale thereof. A copy of this Order may be filed as a Financing Statement with the appropriate state authorities to perfect the security interest of Sellery. 4. CANCELLATION OF PRODUCT. Customer hereby acknowledges and agrees that all Product is custom manufactured for each individual Customer and

once production has commenced there can be no cancellation of orders.

5. TERMINATION. This Order may be terminated by either party at any time: (a) if the other party has breached any material provision contained in this Order and has not cured the breach within thirty (30) days of receipt of written notice of the breach; or (b) if the other party petitions for relief under the Bankruptcy Act, if a receiver is appointed to conduct its business, if its makes an assignment for the benefit of creditors or if an involuntary bankruptcy petition is brought against it and has not been discharged within 30 days of date brought. Seller termination of this Order shall be without prejudice to any other remedies Seller may have. In the event Seller terminates this Order pursuant to Article 4, in addition to any other remedies available to Seller. Customer shall be obligated to pay for all Product shipped and shall accept delivery of and pay for all Product on order as of the date of termination. Customer's obligation to pay all charges accrued as of the termination hereof shall survive such termination. 6. GENERAL PROVISIONS. These terms and conditions shall supersede all terms and conditions printed or implied on any purchase order issued by Customer. No other terms or conditions shall be binding upon Customer or Seller unless said terms and conditions are specifically accepted in writing by Customer and Seller. The failure of Customer or Seller to enforce at any time or for any period of time any of the provisions of this Order shall not constitute a waiver of such provisions. Seller shall not be liable for delays or failure to ship for causes beyond Seller's control, affecting it or its suppliers, including, in part: natural disasters: war: Customer's conduct: labor unrest: riot: governmental interference; unavailability or shortage of labor, materials, fuel or power through customary means at normal rates; failure or destruction of plant or equipment arising from any cause whatsoever; or failure of transportation. If any provision of this Order should be determined to be unenforceable or illegal, it shall be deemed severable from the other provisions. which shall remain valid and enforceable. The validity. performance and construction of this Order shall be governed by the laws of the state of California. Customer agrees to accept venue in, and accept the exclusive jurisdiction of, the courts sitting in the state of California. Any disputes arising out of the performance of this Order shall be settled in accordance with the rules of the American Arbitration Association. The arbitrator may award all costs of dispute resolution, including attorneys' fees, to the prevailing party. Seller agrees that the Product to be furnished under this Order was produced in compliance with all applicable federal requirements with respect to fair labor standards and equal employment opportunity.

7. LIMITED WARRANTY ON PRODUCT. Seller warrants that all Products are warranted to be free of defects in material and workmanship for a period of two (2) years from the date of shipment. This excludes glasses.

Eyewear shipments in the U.S. are covered by a 30-day warranty from date of shipment. Products returned under the above warranty shall be either repaired or replaced at the sole discretion of INFAB, LLC., and is limited to the value of the Products on the date of the original shipment. All replaced parts and Product become the property of Seller. This limited warranty does not include service to repair damage to the Product resulting from customer's negligence, misuse, modification of the Product, or use for purposed other than those intended by its manufacture. Does not cover excessive use or wear and tear from improper storage or customer abuse or failure to follow proper care guidelines. Warranty is voided if any solvent based cleaning solutions are used. Seller's obligations under this limited warranty are conditioned upon Customer's compliance with Seller's reasonable instructions and requests in connection with its performance of such obligations. Service under this limited warranty may be obtained by delivering the Product at Customer's expense (using the original shipping container or its equivalent) during the warranty period to the manufacturer and providing proof of purchase date. Customer shall assume the risk of loss or damage in transit. Seller shall bear the cost for return of repaired or replaced Product under valid warranty claims to Customer by normal commercial means in the continental United States. For further information, refer to the Limited Warranty contained in Seller's standard product literature, if applicable.

NOTE: All aprons must be inspected by the facility before they are put into service. No returns for aprons with alleged pinholes, tears, or cracks will be allowed more than seven (7) days after the aprons have been received by the end user. Protective core material is fragile, easily damaged and therefore great care is needed to protect the integrity of the core material. All aprons are inspected and free of defects when they are shipped. Pinholes, tears, cracks and other damage to the protective core material can be caused by many different factors, but if they are not there when the apron is received and inspected, then they are not the responsibility of the manufacturer and the warranty does not apply.

EXCEPT AS SET FORTH HEREIN AND IN THE LIMITED WARRANTY INCLUDED IN SELLER'S STANDARD PRODUCT LITERATURE, IF APPLICABLE, SELLER HEREBY DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF THE EQUIPMENT, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. OTHER THAN AS STATED ABOVE THE LIMITED WARRANTY SET FORTH IN THIS SECTION, COMBINED WITH THE LIMITED WARRANTIES CONTAINED IN SELLER'S STANDARD PRODUCT LITERATURE, IF APPLICABLE, SHALL CONSTITUTE THE SOLE AND EXCLUSIVE REMEDY OF CUSTOMER AND SHALL ESTABLISH

THE SOLE AND EXCLUSIVE LIABILITY OF SELLER. ITS SUCCESSORS OR ASSIGNS IN CONNECTION WITH THE PRODUCT. IN NO EVENT WILL SELLER BE LIABLE TO CUSTOMER FOR ANY DAMAGES, INCLUDING SPECIAL. INCIDENTAL OR CONSEQUENTIAL DAMAGES. ARISING OUT OF THE USE OR INABILITY TO USE SUCH PRODUCT, EVEN IF SELLER HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR ANY CLAIM BY ANY OTHER PARTY. SOME STATES DO NOT ALLOW THE LIMITATION OR EXCLUSION OF IMPLIED WARRANTIES OR OF LIABILITY FOR SPECIAL OR CONSEQUENTIAL DAMAGES. SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU SOME STATES DO NOT ALLOW THE LIMITATION OF LIABILITY OR THE EXCLUSION OF IMPLIED WARRANTIES. SO THE ABOVE LIMITATIONS AND EXCLUSIONS MAY NOT APPLY AND CUSTOMER MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

8. DAMAGE IN SHIPMENT. The packaging of Infab products meets or exceeds all requirements of the Consolidated Freight Classifications. However, all packages should be checked immediately upon arrival at Customers facility. Should any concealed damage be suspected the Customer shall note it on the Bill of Lading when received. If any concealed damage is discovered, place a hold the shipping container and contact the carrier to schedule an inspection. Any claim must be made directly to the carrier within 24 to 48 hours of merchandise or the claim may not be allowed. Any claim for damage must be made directly with the carrier. Infab cannot file a claim in the name of or for the benefit of the Customer. Failure to follow this exact procedure could invalidate your claim. 9. DELIVERY. All stock items are shipped from our warehouse the same day if P.O. is received before noon. Otherwise it will ship next day. Non-stock or custom items (such as modified or special size aprons and/or monogramming) are shipped within 7 to 10 days of our receipt of your completed order. All special services such as UPS Blue or Red. Fed-Ex. Airborne, etc. can be accommodated at actual cost.

10. MINIMUM ORDERS. There is a \$100.00 minimum per purchase order.


11. EXCHANGES OR RETURNS. Exchanges or Returns must have prior approval by notifying us within 5 days of receipt of the order whether shipped direct to you or drop shipped to your customer. All returns must be sent in their original packaging. A 20% restock charge plus all freight charges including the outbound and return freight will be due and payable. Custom items such as special sized aprons, monogrammed aprons or a product with any modification are not returnable. This includes barriers or racks in any other color except our standard Medical White.




# **View the New Infabcorp.com**

**Effortless Navigation, Enhanced Experience** 








- Buyers Guides
- Resources
- Catalogs
- Blog
- Datasheets
- Certifications

# **Explore the New Custom Apron Configurator**

- Easy-to-use interface
- Real-time customization
- Wide range of colors & patterns
- Add logos & monograms





# **INFAB** 2024

INFAB, LLC info@infabcorp.com infabcorp.com

1040 Avenida Acaso Camarillo, CA 93012 USA Phone: 805.987.5255 Fax: 805.482.8424











REV001 240401