

RADIATION PROTECTION
CATALOG
2026



INFAB

FOR THOSE WHO NEVER STOP

You show up. Every shift. Every patient. Every time.

You deserve protection that works as hard as you do.

TABLE OF CONTENTS

brrr ^o cooling fabric.....	4
UnderShield™.....	5
Revolution Aprons.....	6
Revolution Base Aprons.....	8
Classic Aprons.....	10
Accessories.....	12
Patient Care.....	14
Dental, Drapes, Scatter Protection.....	15
Barriers.....	16
Storage, Curtains, Cleaner.....	17
Smart Track Apron Inspection.....	18
Fit Guide.....	20
Sizing.....	21
About.....	22
Contact.....	23

WHO WEARS INFAB

Neurosurgery | Cardiology | Cath Lab | Radiology | Orthopedics | Oncology | Anesthesiology | Veterinary | Dental



SHOP NOW



THE SCIENCE

Cool under pressure. Literally.

Eight-hour cases. Hot OR lights. The last thing you need is gear that makes you sweat.

brrr° cooling technology uses minerals woven into the fabric to actively pull heat away from your body. Not a coating. Not a gimmick. Real thermal regulation that lasts the life of the garment.



SHOP NOW

brrr° TRIPLE CHILL EFFECT.

HOW IT WORKS



COOLING MINERALS
Instant, continuous cooling effect.
Lowers skin temperature by ~3°F



ACTIVE WICKING
Pulls moisture away from skin for nonstop
comfort during long procedures



RAPID DRYING
Fresh, crisp experience all day.
No more clammy feeling

UNDER APRON



The protection your apron misses.

Traditional aprons can leave gaps. Axillary exposure is real, especially for female surgeons who face 3x higher breast cancer rates.

UnderShield™ provides targeted coverage that works under any apron, with brrr° cooling built in.

KEY FEATURES

Targeted axillary protection

Fully adjustable for a custom fit

Works with any apron brand

brrr° cooling included



SHOP NOW



REVOLUTION

Performance meets style.

Our flagship collection featuring brrr° cooling technology, ergonomic weight distribution, and a variety of colors that look great on you.

Because you shouldn't have to choose between protection and personality.

THE DETAILS

Smooth Revolution Fabric
brrr° Cooling Technology
.25 / .35 / .50mm LE
Lead or Lead-Free Core Protection
Multiple Styles & Colors
XS - 3XL + Custom



[SHOP NOW](#)

VEST & SKIRT STYLES




L103 Lumbar Vest & Skirt



103 Vest & Skirt



703 Full Overlap Lumbar Vest & Skirt



903 Reverse Vest & Skirt

FRONT & FULL PROTECTION



303 Front Protection
403 Front Elastic Back



503 Front Black Belt



803 Black Belt Wrap



603 Full Wrap Elastic Back

18 COLORS



REVOLUTION BASE

Same fit. Your personality.

All the ergonomic design of Revolution in classic nylon.
Over 100 patterns - from planets to floral - and everything
in between.

Professional doesn't have to mean boring.

THE DETAILS

Classic Nylon Fabric
.25 / .35 / .50mm LE
Lead or Lead-Free Core Protection
7 Styles / 100+ Patterns
XS - 3XL + Custom



[SHOP NOW](#)

BASE STYLES



BL103 Lumbar Vest & Skirt Base



B103 Vest & Skirt Base



B203 Full Wrap Base



B303 Front Base



B403 Elastic Back Base



B503 Black Belt Base



B603 Elastic Wrap Base

CUSTOM DESIGNS

Your logo. Your mascot. Your vision. Upload artwork at infabcorp.com or email sales@infabcorp.com

We can match any color and print any design. Hospital logos, team mascots, personal artwork - if you can imagine it, we can make it.

100+ PATTERNS



[VIEW MORE](#)

CLASSIC

Proven. Reliable. Since 1981.

Some things don't need reinventing.
Our Classic collection has protected healthcare
professionals for over four decades.

Every color. Every pattern.
The broadest selection anywhere.

THE DETAILS

Classic Nylon Fabric
.25 / .35 / .50mm LE
Lead or Lead-Free Core Protection
8+ Styles / All Colors
XS - 3XL + Custom



[SHOP NOW](#)

CLASSIC STYLES



VAS Vest & Skirt



RVAS Reverse Vest & Skirt



LSW Lumbar Support Full Wrap



FCW Full Coat Wrap



AWC Weight Reliever Front



EBA Elastic Back Front



LSF Lumbar Support Front



WBT 2" Wide Belt Front

COLORS



PATTERNS



MORE
COLORS &
PATTERNS

ACCESSORIES - COMPLETE PROTECTION

Head to toe. No gaps.

Your apron is the foundation. These are the finishing touches that complete your protection.

AXILLARY



UnderShield™
Targeted axillary protection.
brrr° cooling. XS-3XL.



THYROID



MTC Magnetic
One-hand closure





RADIATION EYEWEAR

FEATURING TRUSTED BRANDS



GUCCI


Ray-Ban

.75mm LE protection. Prescription available. Designer frames from brands you know.

GLOVES

CAPS

MONOGRAMS & PATCHES



Sterile radiation protection.
Lead-free available.

Skull protection. brrr^o cooling.
.25-.50mm LE

Block & script font options.
Adhesive hook & loop patches.

PATIENT CARE

They're in your hands.

From positioning to protection, give every patient, at every age, the care they deserve.


PATIENT PROTECTION

Complete selection of patient protection products. Disposable options also available.



PATIENT POSITIONING

Precision positioning aids for weight-bearing studies, extremities, and full-body imaging. Radiolucent materials won't interfere with image quality.



DENTAL



Adult and pediatric aprons.
Pano capes. Lightweight comfort.



DRAPES



Sterile drapes for procedures.
Multiple configurations.

SCATTER PROTECTION

Scatter Armor shields provide scatter radiation protection for radial and many other types of procedures.



FACILITY

Protect the whole room.

Mobile barriers, curtains, scatter shields.
Room-level protection for real workflows.

MOBILE BARRIERS

2.0mm Pb LE lead glass window
1.5mm Pb LE
Hospital-grade casters with brakes
Standard / Wide / X-Tall / Custom



[REQUEST QUOTE](#)

STORAGE



Mobile Racks



Lockers



CURTAINS



Ceiling and wall mounted.
Custom sizes available.



Wall systems, hangers, racks, lockers. Proper storage extends apron life and keeps your investment protected.

CARE

Always hang - never fold
Wipe after each use
Approved cleaners only
Inspect regularly



TRAVEL BAG

On-the-go Solutions




SMART TRACK

Track. Inspect. Report.

Free garment tracking. Document inspections, set reminders, generate compliance reports. No software. No fees.

FEATURES

- Garment dashboard
- Inspection reminders
- Compliance reports
- Bulk uploads
- Mobile app
- Works with any apron



START FREE

smartrack.infabcorp.com

FROM THE FIELD

"BIG thank you to INFAB because this apron is SERIOUSLY impressive. Lightweight, flexible, and actually comfortable to wear all day, which is saying a LOT if you work in imaging. No stiffness, no dragging you down, just clean protection that moves with you."

- Diego Diaz, RT(R)

"INFAB is the FIGS® of radiation protection. Cool, modern, and built for professionals who care about both performance, comfort, and style. It's not just gear; it's a statement that safety can look good."

- Alvaro Riojas, MD, Pediatrics Cardiologist

"I have been loving the UnderShield! It is a huge help at our trauma hospital where lead in women's sizes can be scarce or in the ED where communal lead is typically extra large and I am standing right next to a mini C-arm during reductions. I get so many compliments..."

- Jessica Schmerler, MD, Orthopedic Surgeon

"The amount of compliments I get is unreal! Thanks to INFAB for keeping me stylish and protected!"

- Mariela Jara, RT(R)

Join 50,000+ healthcare professionals who trust INFAB for their radiation protection.

TRUSTED BY



FIT GUIDE

The right fit changes everything.

Better protection. More comfort. Longer life.
Take five minutes to measure.

HOW TO MEASURE

CHEST (C)

Measure around fullest part

WAIST (W)

Measure at natural waistline

HIPS (H)


Measure widest part of hips

LENGTH (L)

Shoulder to 1" above knee



CHOOSE YOUR SIZE



STANDARD SIZES

Female Sizes

Size	Chest	Waist
XS	30-32"	25-27"
S	33-36"	28-30"
M	37-39"	31-33"
L	40-42"	34-36"
XL	43-45"	37-39"
2XL	46-48"	40-42"

Male Sizes

Size	Chest	Waist
S	34-37"	28-31"
M	38-41"	32-35"
L	42-45"	36-39"
XL	46-49"	40-43"
2XL	50-52"	44-46"

* Models 603 & 903 have unique sizing charts.

TIPS

Measure over light clothing
Between sizes? Go up
Length alterations +/-2" or +/-4"
Custom sizing for perfect fit
Classic aprons for more room.
Revolution aprons for slim fit.

FULL SIZE CHARTS

Detailed male and female size charts with measurements at infabcorp.com/sizing
Or scan QR code for interactive fitting guide.



BOOK A FITTING

NEED HELP?

805-987-5255
sales@infabcorp.com

Real humans. Real help.
We'll get you fitted right.

OUR STORY

American-made. Since 1981.

For over 40 years, we've designed and manufactured radiation protection products in Camarillo, California.

Every apron. Every collar. Made-to-order by people who understand what you do and why protection matters.

**From order to delivery, experience the fastest
turnaround time in the industry.**



40+

Years

100%

Made in USA

1000s

Healthcare Facilities



CERTIFICATIONS

ISO 13485:2016

FDA Registered Medical Device Manufacturer

CE and UKCA conformity for designated products

LET'S TALK

805-987-5255

Monday - Friday 7am - 4pm Pacific

info@infabcorp.com

infabcorp.com

INFAB, LLC
1040 Avenida Acaso
Camarillo, CA 93012

ABOUT US



FOLLOW US



<https://linktr.ee/infabcorp>

Instagram | LinkedIn | Facebook | TikTok | X

SALES INQUIRIES

sales@infabcorp.com
Volume pricing available
Facility outfitting programs

TERMS & CONDITIONS



INFAB

For Those Who Never Stop.



infabcorp.com | 805-987-5255

1040 Avenida Acaso, Camarillo, CA 93012

© 2026 INFAB, LLC. Made in USA