



**UFW Lightning
Full Wrap**



**UVS Lightning
Vest & Skirt**



**UFP Lightning
Front Protection**

These aprons are ready to ship within one week and come with advanced features for enhanced protection and comfort. Available in six vibrant colors, these aprons are designed to meet the radiation protection needs of healthcare professionals.



**Ships within
One Week**



**Unisex
Sizing**



**Six
Colors**



**Smart
Track**



**Abrasion
Resistant**



**Water
Resistant**

Leading the Way in Radiation Protection
Visit our website: www.infabcorp.com | Call us: 805.987.5255

LIGHTNING APRONS BY INFAB

PRODUCT DETAILS, COLORS & SIZING



UVS Lightning Vest & Skirt

The Lightning Vest and Skirt apron provides partial overlap protection and maximum weight distribution between the shoulders and hips, eliminating stress on the upper and lower back. SmartTrack enabled.

- SIZE:** Small - 2xLarge
FABRIC: Nylon (with a double welt pocket)
PROTECTION: .35mm LE, and .50mm LE
CORE: KIARMOR Lead-Free, Lightweight Lead,
COLORS: Black, Blue, Charcoal, Purple, Red, Teal



UFW Lightning Full Wrap

The Lightning Full Wrap offers maximum front and back partial overlap protection. Built-in slits on each side with hook and loop closures provide increased mobility. SmartTrack enabled.

- SIZE:** Small - 2xLarge
FABRIC: Nylon (with a double welt pocket)
PROTECTION: .35mm LE, and .50mm LE
CORE: KIARMOR Lead-Free, Lightweight Lead,
COLORS: Black, Blue, Charcoal, Purple, Red, Teal



UFP Lightning Front Protection

The Lightning Front Protection apron features side closure straps at the waist which relieve weight from the shoulders and adjust for the most comfortable fit. SmartTrack enabled.

- SIZE:** Small - 2xLarge
FABRIC: Nylon (with a double welt pocket)
PROTECTION: .35mm LE, and .50mm LE
CORE: KIARMOR Lead-Free, Lightweight Lead,
COLORS: Black, Blue, Charcoal, Purple, Red, Teal

Add-on: Only Available with Purchase of Lightning Apron



UST-S1 • Thyroid Collar without Binding

Key Features:

A comfortable fit that is fully adjustable with hook and loop closure

* Color options are the same listed above

Colors Available



Sizing Chart

| Size | Chest | Waist | Length |
|------|---------------------|---------------------|-------------|
| S | 34-37" / 86-94 cm | 28-31" / 71-79 cm | 38" / 97 cm |
| M | 38-41" / 97-104 cm | 32-35" / 81-89 cm | 38" / 97 cm |
| L | 42-45" / 107-114 cm | 36-39" / 91-99 cm | 38" / 97 cm |
| XL | 46-49" / 117-125 cm | 40-43" / 102-109 cm | 38" / 97 cm |
| 2XL | 50-52" / 127-132 cm | 44-46" / 112-117 cm | 38" / 97 cm |