



BARRIERS

Medium Window Mobile Barrier

683468



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.

Model Number	Description
683468	Medium Window Mobile Barrier - Lead
Dimensions	74" H x 52½" W X 25" D
Glass Size	12" H x 24" W
Lead Glass	2.0 mm Pb LE
Frame Color	White
Frame	Powder Coated Steel
Floor Curtain Color	Revolution Grey
Wheels	4 swivel casters, 2 with brakes
Crate Size	73" x 30" x 57"
Crate Weight	410 lbs / 185.97 kg



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB