




ACCESSORIES

Zip Top Cassette Covers

BOX OF 100/500



Zip Top



- Zip top
- Packed in handy dispenser boxes
- 100 or 500 individual covers per box
- Available in the five cassette sizes
- Always in stock
- Impervious to fluids
- Available in money saving cases of 500.
- Easy to use
- High quality
- Sanitary
- Disposable
- Top grade polyethylene plastic

Model	Description		Dimensions	Fits Cassette Size
910810-ZT	Cassette Cover Zip Top	Disposable Box of 100	11.25" x 13.25"	8" X 10"
900810-ZT	Cassette Cover Zip Top	Disposable Box of 500	11.25" x 13.25"	8" X 10"
911012-ZT	Cassette Cover Zip Top	Disposable Box of 100	13.375" x 14.75"	10" X 12"
901012-ZT	Cassette Cover Zip Top	Disposable Box of 500	13.375" x 14.75"	10" X 12"
911417-ZT	Cassette Cover Zip Top	Disposable Box of 100	17" x 19.875"	14" X 17"
901417-ZT	Cassette Cover Zip Top	Disposable Box of 500	17" x 19.875"	14" X 17"
911619-ZT	Cassette Cover Zip Top	Disposable Case of 100	19.25" x 22.25"	16" x 19"
901619-ZT	Cassette Cover Zip Top	Disposable Case of 500	19.25" x 22.25"	16" x 19"
912420-ZT	Cassette Cover Zip Top	Disposable Case of 100	20" x 24.625"	24" x 20"
902420-ZT	Cassette Cover Zip Top	Disposable Case of 500	20" x 24.625"	24" x 20"