



GLASSES

GUCCI

RADIATION PROTECTION GLASSES



GG0343-007
Multicolor

The Gucci GG0343 features a rectangular shape. The temples feature signature Gucci color bands and subtle branding.

Temple Length: 150mm



Lens Height: 40mm



Lens Width: 57mm



Bridge Width: 18mm



Frame Size	57mm - 18mm - 150mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Multicolor
Weight	79 G
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

