



GLASSES

GUCCI

RADIATION PROTECTION GLASSES



GG0613-002
Havana Gold




GG0613-003
Pink Transparent

Temple Length: 140mm



Lens Width: 52mm



Bridge Width: 16mm



Feminine Gucci radiation protection eyewear combines luxury design heritage with small bridge styling and 0.75mm lead equivalent protection for elegant clinical safety.

Frame Size	52mm - 16mm - 140mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Havana Gold, Pink Transparent
Weight	TBD
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

