



GLASSES

GUCCI

RADIATION PROTECTION GLASSES



GG1447-001
Black



GG1447-004
Grey Transparent

A best-selling rectangular frame reimaged in a sleek, lightweight acetate. Features refined Gucci logo detailing on the temples and subtle metal accents at the hinges. A timeless design that blends modern minimalism with high-quality craftsmanship.

Lens Width: 57mm

Bridge Width: 17mm



Temple Length: 145mm



| | |
|--------------|---|
| Frame Size | 57mm - 17mm - 145mm Alternative Fit |
| Prescription | Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive |
| Colors | Black, Grey Transparent |
| Weight | 68 G |
| Protection | .75 mm LE |
| Includes | Case, Cloth & Anti-Fog Spray |



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

