



GLASSES

# GUCCI

## RADIATION PROTECTION GLASSES



GG1573-001  
Matte Black

Temple Length: 145mm



Lens Width: 55mm



Bridge Width: 16mm



Lightweight rectangular radiation protection glasses featuring 0.75mm lead equivalency and subtle Gucci branding for sophisticated safety in clinical environments.

<b>Frame Size</b>	55mm - 16mm - 145mm
<b>Prescription</b>	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
<b>Colors</b>	Matte Black
<b>Weight</b>	70 G
<b>Protection</b>	.75 mm LE
<b>Includes</b>	Case, Cloth & Anti-Fog Spray



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

