



GLASSES

GUCCI

RADIATION PROTECTION GLASSES



GG1612-001
Black Crystal

Temple Length: 145mm



Lens Width: 54mm



Bridge Width: 16mm



Distinguished Gucci radiation protection glasses showcase iconic flag colorway temples with adjustable nose pads and 0.75mm lead equivalent lenses, merging luxury with essential safety.

Frame Size	54mm - 16mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Black Crystal
Weight	TBD
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray