



GLASSES

GUCCI

RADIATION PROTECTION GLASSES



GG1737-001
Black

Temple Length: 145mm



Lens Width: 54mm



Bridge Width: 16mm



Bold Gucci rectangular radiation protection eyewear features distinctive emblem detailing and signature flag colorway accents with 0.75mm lead equivalent protection for style-conscious professionals.

Frame Size	54mm - 16mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Black
Weight	67 G
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

