



GLASSES

# Nike 7015 Lead Glasses



**7015-02**  
Matte Black/Space Blue


**Bridge Width:** 17mm



**Temple:** 145mm



**Lens Width:** 55mm



**Lens Height:** 38.1mm



The Nike 7015 INFAB radiation protection glasses combine a classic rectangular acetate frame with 0.75mm lead equivalency lenses, featuring a comfortable single bridge design and durable hinges perfect for healthcare professionals.

<b>FRAME SIZE:</b>	Universal
<b>WEIGHT:</b>	65 g
<b>PRESCRIPTION:</b>	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
<b>COLORS:</b>	Matte Black/Space Blue
<b>PROTECTION:</b>	.75mm LE
<b>INCLUDES:</b>	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

