



GLASSES

Nike 7160 Lead Glasses




7160-01
Clear Forest




7160-02
Crystal Midnight Navy


Bridge Width: 17mm



Temple: 145mm



Lens Width: 55mm



Lens Height: 38.8mm



The Nike 7160 INFAB radiation protection glasses feature a lightweight rectangular design with flexing spring-loaded hinges and texturized rubber temple tips that stay put without pinching. These glasses incorporate 0.75mm lead equivalency lenses while maintaining optimal comfort and a secure fit that complements round and wide face shapes.

FRAME SIZE:	Universal
WEIGHT:	72 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Clear Forest Crystal Midnight Navy
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

