



GLASSES

Nike 7177 Lead Glasses



7177-01
Cacao/Pink Foam Gradient




7177-02
Blue/Peach Gradient

The Nike 7177 radiation protection glasses combine ultra-lightweight construction in high quality acetate with 0.75mm lead equivalency lenses. These glasses deliver optimal stability and comfort during long medical procedures, perfect for professionals who value both style and function.


Bridge Width: 17mm



Temple: 140mm



Lens Width: 53mm



Lens Height: 39.5mm



FRAME SIZE:	Universal
WEIGHT:	67 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Cacao/Pink Foam Gradient Blue/Peach Gradient
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

