



GLASSES

Nike 7405 Lead Glasses



7405-01
Matte Black



7405-02
Matte Obsidian




7405-03
Matte Vintage Green

The Nike 7405 INFAB radiation protection glasses offer a modern angular design with rubberized temple tips for secure fit, featuring 0.75mm lead equivalency lenses in a durable matte-finish frame for essential medical radiation protection.


Bridge Width: 18mm



Temple: 145mm



Lens Width: 56mm



Lens Height: 37mm



FRAME SIZE:	Universal
WEIGHT:	7405-01 - 60 g 7405-02 - 62 g 7405-03 - 60 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Black Matte Obsidian Matte Vintage Green
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case