



GLASSES

Nike Maverick Lead Glasses



MAVERICK-04
Matte Midnight Navy



MAVERICK-03
Matte Cool Grey

The Nike Maverick glasses are pure performance on the inside. The lightweight frame features ventilation and temple arms that provide a soft, firm grip, ensuring you can wear these sunglasses for hours at a time.

Lens Width: 59mm



Lens Height: 45.3mm



Bridge Width: 15mm



Temple Length: 145mm



FRAME SIZE:	Universal
WEIGHT:	78 G
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Midnight Navy, Matte Cool Grey
PROTECTION	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

