



GLASSES

# Maxx 30

**830250** Black**830251** Red**830252** Silver**800070** Blue*Optional*  
Side Protection

- Prescription compatible (for most Rx)
- First protective eyewear to receive the American Ophthalmology Associations seal of acceptance.
- Meets ANSI Z87.1+ standards along with ASTM F803 impact resistance standards.
- Temple bars and nose bridge have soft rubberized padding for extra comfort
- Lightweight and comfortable
- Large lens area for maximum protection
- Bifocal Progressive Rx frame (requires segment height)

Lens Width: 53mm

Bridge Width: 17mm

Temple Length:  
130mm

Lens Height: 35mm

Frame Width: 128mm



Frame Size	Universal
Prescription	Plano Single Vision (-SV) Bifocal (-BF) Bifocal Progressive (-BF-P)
Colors	Black, Red, Silver, Blue
Weight	63G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USAtel: 805.987.5255  
fax: 805.482.8424info@infabcorp.com  
www.infabcorp.com**INFAB**