



GLASSES

# Nike Radeon 3 Lead Glasses



**RADEON3LB-01**  
Shiny Black/Grey



**RADEON3LB-02**  
Clear Blue

**Bridge Width:** 18mm



**Temple:** 145mm



**Lens Width:** 55mm



**Lens Height:** 44.6mm



The Nike Radeon 3 Low Bridge INFAB radiation protection glasses feature a modern square-shaped frame with a single bridge design and 0.75mm lead equivalency lenses, delivering essential protection in a sporty, unisex style. If standard glasses typically slip down your nose or rest on your cheeks, the Low Bridge design features a modified nose pad for a closer fit.

<b>FRAME SIZE:</b>	Universal
<b>WEIGHT:</b>	RADEON3LB-01 - 80 g RADEON3LB-02 - 81 g
<b>PRESCRIPTION:</b>	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
<b>COLORS:</b>	Shiny Black/Grey Clear Blue
<b>PROTECTION:</b>	.75mm LE
<b>INCLUDES:</b>	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

