



GLASSES

# Nike Temper Flair Lead Glasses



**TEMPERFLAIR-01**  
Washed Coral/Mauve



**TEMPERFLAIR-02**  
Vintage Green




**TEMPERFLAIR-03**  
Taupe Grey/Teal

The Nike Temper Flair INFAB radiation protection glasses feature a sophisticated square frame design with 0.75mm lead equivalency lenses for medical professionals. The durable acetate construction and precise single bridge design ensure long-lasting comfort during extended medical procedures.


**Bridge Width:** 18mm



**Temple:** 140mm



**Lens Width:** 53mm



**Lens Height:** 41.4mm



<b>FRAME SIZE:</b>	Universal
<b>WEIGHT:</b>	TEMPERFLAIR-01 - 68 g TEMPERFLAIR-02 - 66 g TEMPERFLAIR-03 - 68 g
<b>PRESCRIPTION:</b>	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
<b>COLORS:</b>	Washed Coral/Mauve Vintage Green Taupe Grey/Teal
<b>PROTECTION:</b>	.75mm LE
<b>INCLUDES:</b>	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

