



GLASSES

Oakley AirDrop



AIRDROP-01
Satin Black

AIRDROP-02
Grey Shadow



Lens Width: 55mm

Bridge Width: 18mm

Lens Height: 33mm



Temple Length: 143mm

Hinge to Hinge: 137mm



The Oakley AirDrop lead glasses feature lightweight O-Matter frame material with a comfortable Three-Point Fit and unobtainium nose pads for secure all-day wear. Fitted with 0.75mm lead equivalent lenses for reliable radiation protection during medical procedures.

FRAME SIZE:	Universal
WEIGHT:	58 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	AIRDROP-01 - Satin Black AIRDROP-02 - Grey Shadow
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

