



GLASSES

Oakley Crosslink Zero



CROSSLINK-01
Matte Black

Lens Width: 55.9mm



Lens Height: 35.4mm



Bridge Width: 16mm



Temple Length: 138mm



The Oakley Crosslink Zero lead glasses blend sporty design elements from the Zero sunglasses with the popular Crosslink eyewear platform. Featuring fixed temples with hidden wire cores for adjustability, Unobtainium nosepads and earsocks for no-slip grip, and 0.75mm lead equivalent lenses for reliable radiation protection during medical procedures.

FRAME SIZE:	Universal
WEIGHT:	61.6 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Black
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

