



GLASSES

Oakley Double Down



DOUBLEDOWN-01
Satin Black



Lens Width: 56mm

Bridge Width: 16mm

Lens Height: 40mm



The Oakley Double Down lead glasses feature a low bridge fit design with lightweight Nano-Matter construction for enhanced comfort. Equipped with 0.75mm lead equivalent lenses for dependable radiation protection during medical and interventional procedures.

FRAME SIZE: Universal

WEIGHT: 71 g

PRESCRIPTION: Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

COLORS: Satin Black

PROTECTION: .75mm LE

INCLUDES: Anti-fog spray, safety strap add-on, protective case

Temple Length: 143mm

Hinge to Hinge: 134mm



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

