



GLASSES

Oakley Holbrook XL



HOLBROOKXL-01
Matte Black



Lens Width: 59mm

Bridge Width: 18mm

Lens Height: 45.5mm



Temple Length: 137mm

Hinge to Hinge: 140.3mm



The Oakley Holbrook XL lead glasses bring the iconic Holbrook design into lightweight O-Matter construction with metal icon accents. Available in standard RX and XL sizes. Fitted with 0.75mm lead equivalent lenses for reliable radiation protection in clinical settings.

FRAME SIZE:	Universal
WEIGHT:	80 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Black
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

