



GLASSES

Oakley Jupiter Squared



JUPITER-01
Matte Black

Lens Width: 56mm



Lens Height: 39.9mm



Bridge Width: 18mm



Temple Length: 135mm



The Oakley Jupiter Squared lead glasses feature bold rectangular styling with switchblade-knife-inspired earstems, metal bands and rivets, and Unobtainium earsocks that grip better with perspiration. O-Matter frame provides all-day comfort. Fitted with 0.75mm lead equivalent lenses for reliable radiation protection during medical procedures.

FRAME SIZE:	Universal
WEIGHT:	74.1 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Black
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

