



GLASSES

# Oakley Leadline



**LEADLINE-01**  
Velvet Black



Lens Width: 54mm



Bridge Width: 17mm



Lens Height: 41.1mm



The Oakley Leadline lead glasses feature a round-edge square lens shape with BiO-Matter frame construction for a comfortable, sustainable design. Equipped with 0.75mm lead equivalent lenses for dependable radiation protection in medical environments.

**FRAME SIZE:** Universal

**WEIGHT:** 64 g

**PRESCRIPTION:** Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive

**COLORS:** Velvet Black

**PROTECTION:** .75mm LE

**INCLUDES:** Anti-fog spray, safety strap add-on, protective case

Temple Length: 139mm



Hinge to Hinge: 133mm



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

