



GLASSES

Oakley Turbine



TURBINE-01
Matte Black



Lens Width: 65mm



Bridge Width: 17mm



Lens Height: 43mm



Temple Length: 132mm



Hinge to Hinge:
137mm



The Oakley Turbine lead glasses combine a bold wraparound design with lightweight O-Matter frame material and Unobtanium nose pads for a secure fit. Fitted with 0.75mm lead equivalent lenses for reliable radiation protection during interventional procedures.

FRAME SIZE:	Universal
WEIGHT:	86 g
PRESCRIPTION:	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
COLORS:	Matte Black
PROTECTION:	.75mm LE
INCLUDES:	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

