




GLASSES

Ray-Ban 4378



4378-01
Havana

Temple Length: 145mm



Lens Height: 45.2mm



Lens Width: 52mm



Bridge Width: 16mm



Frame Size	52mm - 16mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Havana
Weight	67 g
Protection	.75mm LE
Includes	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

