




GLASSES

Ray-Ban 5428



5428-01
Havana on Light Blue

Temple Length: 145mm



Lens Height: 36.6mm



Lens Width: 53mm



Bridge Width: 17mm



Frame Size	53mm - 17mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Havana on Light Blue
Weight	64 g
Protection	.75mm LE
Includes	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

