



GLASSES

Ray-Ban 8908




8908-01
Transparent Blue

8908-02
Transparent Grey



Temple Length: 145mm



Lens Height: 38.9mm



Lens Width: 53mm



Bridge Width: 18mm



Frame Size	53mm - 18mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Transparent Blue, Transparent Grey
Weight	64 g
Protection	.75mm LE
Includes	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

