



GLASSES

Ray-Ban New Wayfarer



WAYFARER-01
Rubber Black

WAYFARER-02
Blue



Temple Length: 145mm

Lens Height: 40.8mm



Lens Width: 52mm

Bridge Width: 18mm



Frame Size	52mm - 18mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Rubber Black, Blue
Weight	71 g
Protection	.75mm LE
Includes	Anti-fog spray, safety strap add-on, protective case



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

