



# Adult Oblique Finger Block - Coated

## INF-NCAA



Coated Sponges feature a patented medical grade, fluid-resistant coated surface that is easy to clean with any medical grade disinfectant wipe or spray. Coated surface does not easily peel off or tear like other coated sponges on the market. They are radiolucent and have minimal shadowing with low attenuation. These unique heavy-duty sponges are designed and manufactured with a patented flared angle that significantly reduces artifacts during x-ray imaging. The flared sides increase the overall size of a standard non-flared sponge but does not affect imaging or positioning. A superior design for superior images.

<b>Model</b>	INF-NCAA
<b>Description</b>	Adult Oblique Finger Block - Coated
<b>Color</b>	Black
<b>Dimensions</b>	10.5" x 6" x 3"



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

**INFAB**