


PATIENT

Rectangle Block Positioner - Closed Cell

INF-CCBN



Closed Cell Sponges are made from a lightweight, radiolucent, non-porous foam. These positioning sponges are entirely fluid resistant, easy to clean and extremely durable. Featured in a non-coated, premium charcoal grey color.

- 100% Radiolucent
- Artifact free imaging
- Cleanable, fluid resistant surface
- Lightweight and comfortable

Model	INF-CCBN
Description	Closed Cell Rectangle
Color	Charcoal Grey
Dimensions	(LxWxH) 18" x 14" x 3"