



PATIENT

Accessory Rails - Siemens

INF-30465



The INFAB Accessory Rails are available in various lengths and designed to securely clamp to the tabletop and support various surgical accessories. The rails enhance operating room efficiency by enabling quick conversions, thanks to their easy setup and teardown. They can be easily cleaned, making them a practical and reliable addition to any operating room.

- 40" detachable accessory rail for Siemens
- Size: 9.5mm x 28.6mm - set of 2
- Designed to securely clamp to the tabletop and support various accessories useful during surgery such as Retractors, IV Poles, or Armboards.
- Increases operating room efficiency by allowing for quick conversion of operating rooms.
- Quick and easy setup and teardown
- Easily cleaned

Model	INF-30465
Description	Accessory Rails
Modalities	Cath Lab, Hybrid OR, Interventional, Interventional X-Ray, Operating Room, Vascular
Substrate	Aluminum / Stainless Steel
Dimensions	9.5mm x 28.6mm
Finish	ScanCoat®
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB