



Philips AD7 Mayfield Carbon Fiber Extender - Neuro Table

INF-29632



INFAB's AD7 Neuro Carbon Fiber Extender for the Philips AD7 Neuro Table extends the length of the table. It allows for the attachment of Mayfield (Infinity XR2) and Doro (Lucent) cranial fixation devices. This table extender is the perfect addition to your Hybrid OR, with its quick setup and takedown procedures allowing your OR to quickly adapt to changing patient needs.

- Carbon fiber extender for the Philips AD7 Neuro table
- Includes two 15" clamp-on rails
- Allows for attachment of Mayfield (Infinity XR2) and Doro (Lucent) cranial fixation devices
- Allows for full c-arm rotation
- Removable 15" American rails
- Quick and easy setup and takedown
- Radiolucent
- Easily cleaned

Model	INF-29632
Description	Philips AD7 Mayfield Carbon Fiber Extender - Neuro Table
Modalities	Cath Lab, Hybrid OR, Interventional, Interventional X-Ray, Operating Room, Vascular
Substrate	Carbon Fiber
Dimensions	52" x 21.5" x 0.22" Shipping: 58" x 24" x 7"
Weight	25 lbs
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB