




PATIENT

Aquilion Table Pad Cover

INF-9121



These Table Pad covers to protect your Aquilion Table pads from contamination. These covers are vinyl "sleeves," which can easily be slid onto and removed from the table pad. These covers have cutouts in the bottom which preserve the ability of the table pad to Velcro to the Aquilion table. The covers are also disposable.

- Covers Aquilion CT table pads 9106F, 9107, 9108, and 9116
- Provides a barrier for the table pad from contamination
- Frosted vinyl cover "sleeve" with flap
- Design allows for quick and easy installation and disposal
- Access holes on the bottom side of the cover allow adherence to tabletop Velcro
- Single-Use
- Non-Sterile
- Latex Free

Model	INF-9121
Description	Aquilion Table Pad Cover
Modalities	CT Scan, General Positioning, Nuclear Medicine, Ultrasound, X-Ray
Substrate	CT Scan
Packaging	Box of 4
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB