



Canon Aquilion Table Pad w/Fluid Dam

INF-9108



The Canon Aquilion Table Pad is a closed-cell foam pad designed for the Canon Aquilion CT machine. This table pad is also used in tandem with disposable vinyl covers (INF-9121) to increase cleanability and reduce contamination, maximizing the life of the product. This table pad is a perfect addition to your CT machine, maximizing compatibility and patient comfort.

- Closed cell foam table pad for use with the Aquilion CT Tables.
- Ensures rigidity, durability, and fluid resistance.
- Separate vinyl covers (INF-9121) are available to increase cleanability and reduce contamination.
- Closed cell foam is rigid, fluid-resistant, and durable.
- Latex Free
- Reusable
- Radiolucent

Model	INF-9108
Description	Canon Aquilion Table Pad w/ Fluid Dam
Modalities	CT Scan
Substrate	Closed cell foam
Dimensions	1" H x 18.6" W x 96.30" L
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

