



PATIENT

Circular Disc with ScanCoat®

INF-112-SCB



The INFAB Circular Disc is the perfect way to raise a patient's head off the table to ensure proper spinal positioning. Constructed with an open cell foam substrate, it ensures optimal support and comfort for the patient. The disc can is coated with ScanCoat®, a sprayed coating that significantly improves cleanability and durability while maintaining complete radiolucency. The disc is latex-free and reusable, offering a cost-effective and sustainable solution for medical facilities.

- Sprayed with ScanCoat® coating which improves the cleanability and durability of the product, while ensuring complete radiolucency. Allowing you to get a clear scan with fewer rescans needed.
- Latex Free
- Reusable

Model	INF-112-SCB
Description	Circular Disc with ScanCoat®
Modalities	CT Scan, General Positioning, MRI, Nuclear Medicine, Ultrasound
Substrate	Open cell foam sprayed with ScanCoat®
Dimensions	1.5"H x 7" Diameter
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

