



PATIENT

Forehead Strap with Adult Pad

INF-9007



INFAB's Forehead Strap includes a foam pad which provides extra comfort for the patient while maintaining proper restraint.

*Please note: Covered items (like this one) are not recommended for low KVP techniques as the edges/seams may show.

- Patient weighted arm support for Infinix tables.
- Support slides under the table pad to be held in place by the patient's weight while positioning their arms.
- Quick and easy setup and teardown
- Easily cleaned
- Reusable
- Radiolucent

Model	INF-9007
Description	Forehead Strap with Adult Pad
Modalities	CT Scan
Substrate	Open cell foam
Finish	Nylon/ScanCoat
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB