



Medvac Immobilizer Hand Pump

INF-27625



The INFAB MedVac Immobilizer Hand Pump is a manually operated hand pump designed to remove air from both MedVac Immobilizers (parts VMR33PUBB & VMR38X01). Although an electrical pump can be used with the MedVac Immobilizers, the Hand Pump is a more portable and reliable option, as it can be used even in areas without wall sockets. When the immobilizers are vacuumed, they become stiff, which allows them to properly immobilize the patient during their MRI scan. The pump is latex-free and reusable, making it a cost-effective and sustainable option for your diagnostic imaging suite. These features together make the MedVac Immobilizer Hand Pump the perfect product to pair with your MedVac Immobilizer.

- Manually operated hand pump designed to remove air from the MedVac Immobilizers (parts VMR33PUBB & VMR38X01).
- Hand pump allows for MedVac Immobilizers to be used in areas without electricity, as the pump serves as a replacement for an electrical air pump.
- When MedVac Immobilizers are vacuumed, they become stiff, allowing for proper immobilization of the patient.
- Latex-Free
- Reusable

Model	INF-27625
Description	Medvac Immobilizer Hand Pump
Modalities	MRI
Substrate	Plastic
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

