



# MRI Comfort Pillow

## INF-463



The INFAB MRI Comfort Pillow is specifically designed to keep patients comfortable during procedures that require the arms to be overhead. Featuring a high-density foam substrate, the pillow ensures proper support and comfort, while its durable nylon cover enhances durability, cleanability, and contamination control. For added comfort, the pillow includes two reinforced hand straps for the patient to hold during the procedure. Additionally, it features supply pouches on the back, allowing for quick and easy access to equipment. The pillow is latex-free, radiolucent for clear imaging, and reusable, making it a practical and efficient solution for enhancing patient comfort during MRI procedures.

- Pillow designed to keep patients comfortable during procedures where the arms must be overhead
- High-density foam substrate ensures proper support and comfort, durable nylon cover increases durability, cleanability, and contamination control.
- The pillow includes two reinforced hand straps for the patient to hold during the procedure for added comfort.
- This item also features supply pouches on the back that allow for quick and easy access to equipment.
- Latex Free
- Radiolucent
- Reusable

<b>Model</b>	<b>INF-463</b>
<b>Description</b>	Leg Rest with Knee Cutouts
<b>Modalities</b>	CT Scan, General Positioning, MRI, Nuclear Medicine
<b>Dimensions</b>	8.75" H x 16" W x 25" L
<b>Substrate</b>	High-density foam with durable nylon cover
<b>Sterile</b>	N
<b>Latex Free</b>	Y
<b>Single Use</b>	N



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

