




Philips AD7 Cath Table - Foot End Extension

INF-33343



INFAB's AD7 Carbon Fiber Cath Table Foot-End Extender for the Philips AD7 Cath Table extends the length of the Philips AD7 table by 24". The extender can hold up to 50 lbs, easily supporting guide wires, catheters, and instruments, allowing you to reduce room congestion caused by free-standing carts and tables. The extender is also quick and easy to attach, and can be done by a single person, allowing your Hybrid OR to rapidly adapt to changing patient needs.

- Extends the Foot-End of the Philips AD7 Cath table by 24"
- Holds a maximum weight of 50 lbs.
- Built-in handle for quick one-person installation
- Easily supports guide wires, catheters, and instruments while reducing room congestion caused by additional free-standing tables or carts
- Setup and takedown are quick and easy
- Latex free
- Easily cleaned and disinfected

Model	INF-33343
Description	Philips AD7 Cath Table - Foot End Extension
Modalities	Cath Lab, Hybrid OR, Interventional, Interventional X-Ray, Operating Room, Vascular
Dimensions	25" x 24" x 0.28"
Weight	14 lbs
Sterile	N
Latex Free	Y
Single Use	N

