



PATIENT

Philips AD7 Neuro Table Pad

INF-33535



The INFAB Table Pad for the Philips AD7 Neuro Table is a full-length, premium foam pad specifically designed for the Philips AD7 Neuro tabletop. Featuring sonic-sealed seams, it offers superior protection against fluid ingress. Incorporating Tempur-pedic™ foam technology, this pad ensures optimal patient comfort and has been developed in conjunction with healthcare professionals. The pad is radiolucent, making it suitable for imaging procedures, and is both latex-free and reusable, ensuring safety and sustainability. Additionally, it is easily cleaned, providing a practical solution for maintaining hygiene in medical settings.

- Full-length, premium foam table pad designed for the Philips AD7 Neuro tabletop
- Head & shoulder shaped design
- Sonic-sealed seams protect against fluid ingress
- 3" foam stacked pad with comfort top layer
- Antifungal/Antimicrobial fluid resistant cover
- PVC and Latex Free
- Radiolucent
- Easily cleaned

Model	INF-33535
Description	Philips AD7 Neuro Table Pad
Modalities	Cath Lab, Hybrid OR, Interventional, Interventional X-Ray, Operating Room, Vascular
Dimensions	3"x22"x125"
Weight	17 lbs
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

