


# Premium Heavy Duty Table Pad

## INF-520



The INFAB Premium Heavy Duty Table Pad is designed for high-use medical environments, providing exceptional comfort and durability. This pad features an open-cell foam substrate with a heavy-duty vinyl cover, enhancing cleanability and contamination control. The open cell foam offers superior comfort and support, ideal for patients lying down for extended periods. Convenient features include grommets for wall storage and a centered crosshair graphic for easy patient positioning. Latex-free and reusable, this pads offers a sustainable, cost-effective solution for healthcare facilities. With a robust design and thoughtful features, it is an excellent choice for enhancing patient care and operational efficiency.

- Open cell foam pad with heavy-duty vinyl cover.
- Heavy-duty vinyl cover improves cleanability and contamination control.
- Covered items like this is not recommended for low KVP techniques as the edges/seams may show.
- Open cell foam substrate improves patient comfort while providing support, making this table pad the perfect solution for situations in which patients must remain lying down for extended periods.
- Grommets are located at the top for wall storage.
- Centered crosshair graphic helps to easily position the patient.
- Latex-Free
- Reusable

<b>Model</b>	INF-520
<b>Description</b>	Premium Heavy Duty Table Pad
<b>Modalities</b>	CT Scan, General Positioning, Nuclear Medicine, Ultrasound, X-Ray
<b>Substrate</b>	Open cell foam
<b>Dimensions</b>	1" H x 24" W x 72" L
<b>Sterile</b>	N
<b>Latex Free</b>	Y
<b>Single Use</b>	N



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel : 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

**INFAB**