



Restraint Belt

INF-168



The INFAB Restraint Belt is a nylon strap with hook and loop fastener lining its entire length. This strap can be used to secure the patient to a table or wheelchair and can be tightened to ensure proper immobilization, security, and comfort. The strap can also be washed when needed. Made from a durable nylon substrate, the strap ensures both longevity and patient comfort. Additionally, the strap is latex-free, radiolucent for clear imaging, fluid-resistant, and reusable, making it a practical and reliable solution for enhancing the quality and efficiency of cervical vertebrae imaging procedures. This strap improves patient safety during operations or transfers, as it can secure the patient to the chair or wheelchair, ensuring that they will be always held in place securely.

- Nylon strap with hook and loop fastener through its entire length.
- Can be used to secure the patient to a table or wheelchair.
- Can be washed.
- Latex Free
- Radiolucent
- Fluid-resistant
- Reusable

Model	INF-168
Description	Restraint Belt
Modalities	General Positioning
Substrate	Nylon
Dimensions	5" W x 72" L
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB