



Steady-Ped™ CT Foot Positioner

INF-33532



The innovative INFAB Steady-Ped™ foot immobilizer helps produce flawless images on the first CT scan. The specialized design secures the patient's feet at a 90° angle with a customizable hook & loop strap. Additional hook & loop is provided for securing the Steady-Ped to the end of the CT table when necessary. The unique curved bottom allows wide compatibility with CT machines, surpassing competitors in precision and performance.

- Secures patient's feet in an upright position (90°) to allow CT imaging of lower extremities and run-offs
- Widely compatible with GE, Canon, and other CT Systems
- The Steady-Ped™ fits in CT machines easily due to its rounded bottom
- The unique curved bottom ensures a perfect fit on CT Machines
- Portable, lightweight, low profile, and reusable
- Easy to use, clean and store
- Radiolucent

Model	INF-33532
Description	Steady-Ped™ CT Foot Positioner
Modalities	CT Scan
Substrate	ABS Rigid Plastic
Dimensions	16" x 11.5" x 11.5"
Weight	3 lbs
Sterile	N
Latex Free	Y
Single Use	N