



# Wedge - 18.75°

INF-105



INFAB's 18.75° Wedge is expertly crafted for optimal general positioning, providing both robust support and patient comfort through its open cell foam substrate. This latex-free, reusable wedge is designed to meet the highest standards of medical safety and efficiency and comes in a set of two.

- 18.75° wedges for general positioning.
- Open cell foam substrate ensures support while providing comfort for the patient.
- Latex Free
- Reusable

<b>Model</b>	INF-105
<b>Description</b>	Wedge - 18.75°
<b>Modalities</b>	CT Scan, General Positioning, MRI, Nuclear Medicine, Ultrasound
<b>Substrate</b>	Open cell foam
<b>Dimensions</b>	3.75"H x 9"W x 8"L
<b>Sterile</b>	N
<b>Latex Free</b>	Y
<b>Single Use</b>	N



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

**INFAB**