


White EZ Strap

INF-135



The INFAB EZ Straps are versatile positioning aids designed to enhance patient stability during medical procedures. Featuring a reliable hook-and-loop fastener system, these straps can be easily cut to the desired length, allowing for precise and customized patient positioning. This adaptability improves their versatility, making them suitable for a wide range of medical applications. They are also washable, ensuring they maintain high hygiene standards and can be reused multiple times. The straps are latex-free, making them safe for patients with latex sensitivities, and are radiolucent, ensuring they do not interfere with imaging procedures. Additionally, their durable construction ensures long-lasting performance, making the INFAB EZ Straps a practical, cost-effective solution for medical facilities seeking to enhance patient care and operational efficiency.

- Hook and loop fastener which may be cut to length to aid patient position
- Can be washed to ensure reusability.
- Ability to cut to desired length improves versatility.
- Latex Free
- Radiolucent
- Reusable

Model	INF-135
Description	White EZ Strap
Modalities	CT Scan, General Positioning, MRI, Nuclear Medicine
Substrate	Hook & Loop
Dimensions	1.5" W x 30' L
Sterile	N
Latex Free	Y
Single Use	N



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB