



Gonad Shields

PGS



- Set of three or available individually: Infant, Youth and Adult
- 3/4" buckle strap
- Easy fastening around the waist
- .50mm LE
- Available in our Classic Nylon Colors

| Model | Description |
|--------|---|
| PGS-S1 | 3" x 5" Infant Size Gonad Shield |
| PGS-M1 | 6" x 7" Youth Size Gonad Shield |
| PGS-L1 | 7" x 9" Adult Size Gonad Shield |
| PGS-X1 | Set of Gonad Shield (Infant, Youth and Adult sizes included) |



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

INFAB