




# Gonadial Diaper

## DPR



Full front and rear protection of a patients genital region in an easy to use and easy to clean diaper style garment. Diapers are available in various sizes that offer a full 0.50mm Pb equivalent of radiation protection. One inch wide buckle closure straps are supplied in addition to hook and loop closure for fast, easy and secure attachment to the patient.

Model	Description
DPR-S	Lead Lined Diaper (Small)
DPR-M	Lead Lined Diaper (Medium)
DPR-L	Lead Lined Diaper (Large)
DPR-X	Lead Lined Diaper (X-Large)

	A	B	C	D
Small	18"	3.5"	23.5"	26"
Medium	18.25"	3.75"	24.5"	30.5"
Large	21"	3.75"	28"	31.5"
X-Large	23"	3.75"	28.25"	30.75"