



Half Apron with Hook and Loop Closure

HSV



- Easy on/off convenience
- Comfortable, durable, patient friendly design
- Traditional hook and loop Velcro type closures
- Wide two inch webbing straps
- Positive and secure closure
- Maximum patient protection
- Three sizes or custom made on request
- Available in your choice of color
- .5 mm Pb protection comes standard

Model Number	Description
HSV-S1	12" x 24" Half Apron with Velcro Closure
HSV-M1	18" x 24" Half Apron with Velcro Closure
HSV-L1	24" x 24" Half Apron with Velcro Closure